


Medical Device Alert

Ref: MDA/2014/040 Issued: 09 October 2014 at 14:00

Device
<p>LIFEPAK[®] 1000 defibrillator.</p> <p>Manufactured by Physio-Control.</p> <p>All serial numbers.</p> 

Problem	Action
<p>Risk of failure to deliver a shock. This is due to unexpected shut down of the defibrillator caused by a battery with very low charge.</p> <p>Confusing instructions for use mean that batteries aren't always replaced when they have a low or very low charge.</p> <p>See the Field Safety Notice for further details.</p>	<ul style="list-style-type: none"> • Identify all LIFEPAK 1000 defibrillators in your possession. • Ensure users read the important information in the Field Safety Notice and follow the instructions on how to check if a new battery is required. • Make sure that a spare, fully charged battery is always available.
Action by	
<p>All staff responsible for the use, maintenance and purchase of these devices.</p>	
CAS deadlines	Contact
<p>Action underway: 23 October 2014</p> <p>Action complete: 06 November 2014</p> <p>Note: These deadlines are for operators to be aware of instructions and to check batteries.</p>	<p>Manufacturer</p> <p>Olivier Anki Physio-Control Operations Tel: +33 695 475 522 Email: olivier.anki@physio-control.com</p>

Distribution

This MDA has been sent to:

- Care Quality Commission (CQC) (headquarters) for information
- Clinical commissioning groups (CCGs)
- HSC trusts in Northern Ireland (chief executives)
- Local authorities in Scotland (equipment co-ordinators)
- NHS boards and trusts in Wales (chief executives)
- NHS boards in Scotland (equipment co-ordinators)
- NHS England area teams for information
- NHS trusts in England (chief executives)
- Special health authorities for information

Onward distribution

Please bring this notice to the attention of relevant employees in your establishment. Below is a suggested list of recipients.

Trusts

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- A&E consultants
- A&E departments
- A&E directors
- A&E nurses
- Adult intensive care units
- Ambulance services directors
- Ambulance staff
- Biomedical engineering staff
- Cardiology departments
- Cardiology nurses
- Cardiology, directors of
- Cardiothoracic departments
- Clinical governance leads
- EBME departments
- Health and safety managers
- Intensive care units
- Intensive care, directors of
- Maintenance staff
- Medical directors
- Medical physics departments
- Nursing executive directors
- Resuscitation officers and trainers
- Risk managers

Independent distribution

Establishments registered with the Care Quality Commission (CQC) (England only)

This alert should be read by:

- Hospitals in the independent sector
- Independent treatment centres

Please note: CQC and OFSTED do not distribute these alerts. Independent healthcare providers and social care providers can sign up to receive MDAs directly from the Department of Health's Central Alerting System (CAS) by sending an email to: safetyalerts@dh.gsi.gov.uk and requesting this facility.

Contacts

Manufacturer

Olivier Anki
Physio-Control Operations Netherlands B.V.
Galjoenweg 68, 6222 NV
Maastricht
The Netherlands
Tel: +33 695 475 522
Email: olivier.anki@physio-control.com

England

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2014/040** or **2014/007/010/081/006**

Technical aspects

Paul Sandhu or Yasser Zayni
Medicines & Healthcare Products Regulatory Agency
Floor 4
151 Buckingham Palace Road
London SW1W 9SZ

Tel: 020 3080 7266 / 7238

Fax: 020 8754 3965

Email: paul.sandhu@mhra.gsi.gov.uk
yasser.zayni@mhra.gsi.gov.uk

Clinical aspects

Mark Grumbridge
Medicines & Healthcare Products Regulatory Agency
Floor 4
151 Buckingham Palace Road
London SW1W 9SZ

Tel: 020 3080 7128

Fax: 020 8754 3965

Email: mark.grumbridge@mhra.gsi.gov.uk

How to report adverse incidents

Please report via our website <http://www.mhra.gov.uk>

Further information about **CAS** can be found at <https://www.cas.dh.gov.uk/Home.aspx>

Northern Ireland

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.

Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre
Health Estates Investment Group
Room 17
Annex 6
Castle Buildings
Stormont Estate,
Dundonald
BT4 3SQ

Tel: 02890 523 704

Fax: 02890 523 900

Email: NIAIC@dhsspsni.gov.uk

<http://www.dhsspsni.gov.uk/index/hea/niaic.htm>

How to report adverse incidents in Northern Ireland

Please report directly to NIAIC, further information can be found on our website <http://www.dhsspsni.gov.uk/niaic>

Further information about **SABS** can be found at <http://sabs.dhsspsni.gov.uk/>

Scotland

All requests regarding return, replacement or modification of the devices mentioned in this alert should be directed to the relevant supplier or manufacturer.

Other enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre
NHS National Services Scotland
Gyle Square
1 South Gyle Crescent
Edinburgh
EH12 9EB

Tel: 0131 275 7575

Fax: 0131 314 0722

Email: nss.irc@nhs.net

<http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-irc/>

Wales

Enquiries in Wales should be addressed to:

Improving Patient Safety Team
Medical Directorate
Welsh Government
Cathays Park
Cardiff
CF10 3NQ

Email: improvingpatientsafety@wales.gsi.gov.uk

MHRA is an executive agency of the Department of Health

© Crown Copyright 2014

Addressees may take copies for distribution within their own organisations