

# Medical Device Alert

Ref: MDA/2013/070 Issued: 28 August 2013 at 15:30

## Device

Insulin infusion sets and reservoirs used with Paradigm ambulatory insulin pumps (manufactured by Medtronic):

- 1) Paradigm infusion sets (all MMT models, all lot numbers).
- 2) Paradigm insulin reservoirs (MMT models - 326A and 332A, specific lots).

See below for details of model and lot numbers for both products.

Problem	Action
<p>Risk of under- or over-delivery of insulin, which may lead to hypoglycaemia or hyperglycaemia, with a consequent loss of consciousness or death, in extreme cases.</p> <p>The root causes are:</p> <ol style="list-style-type: none"> <li>1) Spillage from the insulin reservoir may block connector vents, preventing correct priming of the pump. Medtronic have issued advice on the correct technique for filling the reservoir.</li> <li>2) The insulin reservoir may leak, due to a manufacturing defect, which may result in under delivery of insulin and failure of the pump's occlusion alarm. Medtronic have initiated a <b>recall</b> of affected reservoirs.</li> </ol>	<p>Identify affected infusion sets and reservoirs.</p> <p>Ensure that all relevant staff and patients receive the appropriate copy of Medtronic's <a href="#">Field Safety Notices</a> and are aware of the recommendations. In particular:</p> <ul style="list-style-type: none"> <li>• Ensure staff and patients follow the manufacturer's recommendations to avoid spillage from the insulin reservoir.</li> <li>• Discard unused, affected reservoirs immediately and request replacements.</li> <li>• Return Field safety Action sheet to manufacturer.</li> </ul>
Action by	
<p>Diabetes departments. Those involved in the supply and use of these devices.</p>	
CAS deadlines	Contact
<p>Action underway: 18 September 2013</p> <p>Action complete: 02 October 2013</p> <p><b>Note: These deadlines are for staff and patients to follow the recommendations in the FSNs and to ensure that affected reservoirs are removed from use.</b></p>	<p><b>Manufacturer</b> Lisa Stanley Medtronic Ltd Tel: 01923 212 213 Email: <a href="mailto:Lisa.stanley@medtronic.com">Lisa.stanley@medtronic.com</a></p>

## Device

Models of Paradigm infusion sets affected (all lots):

MMT-317	MMT-318	MMT-324	MMT-325
MMT-312S	MMT-312L	MMT-386	MMT-387
MMT-394	MMT-396	MMT-397	MMT-398
MMT-399	MMT-377	MMT-378	MMT-381
MMT-382	MMT-383	MMT-384	MMT-368
MMT-862	MMT-864	MMT-866	MMT-874
MMT-876	MMT-884	MMT-886	MMT-921
MMT-923	MMT-925	MMT-941	MMT-943
MMT-945	MMT-961	MMT-963	MMT-965
MMT-975			

Lot numbers of recalled MMT-326A and MMT-332A reservoirs:

H8416432	H8461538	H8489386	H8512826
H8420977	H8463297	H8491921	H8515317
H8422490	H8464121	H8492449	H8517079
H8424676	H8467888	H8494645	H8521052
H8437486	H8469703	H8496561	H8539013
H8441420	H8471745	H8500423	H8541843
H8442973	H8473106	H8500472	H8584244
H8451531	H8473271	H8503372	H8627745
H8452933	H8476270	H8503728	H8603292
H8455959	H8478398	H8509305	H8604958
H8457716	H8485398	H8510440	H8635301
H8459557	H8486688	H8512566	

## Distribution

This MDA has been sent to:

- Care Quality Commission (CQC) (headquarters) for information
- Clinical commissioning groups (CCGs)
- HSC trusts in Northern Ireland (chief executives)
- Local authorities in Scotland (equipment co-ordinators)
- NHS boards and trusts in Wales (chief executives)
- NHS boards in Scotland (equipment co-ordinators)
- NHS England area teams
- NHS trusts in England (chief executives)
- Social services in England (directors)

### Onward distribution

Please bring this notice to the attention of relevant employees in your establishment. Below is a suggested list of recipients.

#### Trusts

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- All wards
- Clinical governance leads
- Community children's nurses
- Community diabetes specialist nurses
- Community hospitals
- Diabetes clinics/outpatients
- Diabetes nurse specialists
- Diabetes, directors of
- Diabetologists

- EBME departments
- Equipment stores
- Hospital pharmacies
- Medical directors
- Medical libraries
- Nursing executive directors
- Paediatric nurse specialists
- Paediatricians
- Risk managers
- Supplies managers

### **NHS England area teams**

CAS liaison officers for onward distribution to all relevant staff including:

- Community pharmacists
- District nurses
- General practitioners

### **Social services**

Liaison officers for onward distribution to all relevant staff including:

- Care at home staff
- Care management team managers
- Community care staff

### **Independent distribution**

#### **Establishments registered with the Care Quality Commission (CQC) (England only)**

This alert should be read by:

- Clinics
- Hospitals in the independent sector
- Independent treatment centres
- Private medical practitioners

Please note: CQC and OFSTED do not distribute these alerts. Independent healthcare providers and social care providers can sign up to receive MDAs directly from the Department of Health's Central Alerting System (CAS) by sending an email to: [safetyalerts@dh.gsi.gov.uk](mailto:safetyalerts@dh.gsi.gov.uk) and requesting this facility.

## **Contacts**

### **Manufacturer**

Lisa Stanley  
Medtronic Ltd  
Building 9, Croxley Green Business Park  
Watford, WD18 8WW

Tel: 01923 212 213

Fax: 01923 241 004

Email: [Lisa.stanley@medtronic.com](mailto:Lisa.stanley@medtronic.com)

## **England**

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2013/070** or **2013/006/013/291/041**

### **Technical aspects**

Enitan Taiwo and Claire Dunne  
Medicines & Healthcare Products Regulatory Agency  
Floor 4  
151 Buckingham Palace Road  
London SW1W 9SZ

Tel: 020 3080 7122 / 7162

Fax: 020 8754 3965

Email: [enitan.taiwo@mhra.gsi.gov.uk](mailto:enitan.taiwo@mhra.gsi.gov.uk)  
[claire.dunne@mhra.gsi.gov.uk](mailto:claire.dunne@mhra.gsi.gov.uk)

**Clinical aspects**

Mark Grumbridge  
Medicines & Healthcare Products Regulatory Agency  
Floor 4  
151 Buckingham Palace Road  
London SW1W 9SZ  
Tel: 020 3080 7128  
Fax: 020 8754 3965  
Email: [mark.grumbridge@mhra.gsi.gov.uk](mailto:mark.grumbridge@mhra.gsi.gov.uk)

**How to report adverse incidents**

Please report via our website <http://www.mhra.gov.uk>  
Further information about **CAS** can be found at <https://www.cas.dh.gov.uk/Home.aspx>

**Northern Ireland**

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.  
Enquiries and adverse incident reports in Northern Ireland should be addressed to:  
Northern Ireland Adverse Incident Centre  
Health Estates Investment Group, Room 17, Annex 6, Castle Buildings, Stormont Estate,  
Dundonald BT4 3SQ  
Tel: 02890 523 704 Fax: 02890 523 900 Email: [NIAIC@dhsspsni.gov.uk](mailto:NIAIC@dhsspsni.gov.uk)  
<http://www.dhsspsni.gov.uk/index/hea/niaic.htm>

**How to report adverse incidents in Northern Ireland**

Please report directly to NIAIC, further information can be found on our website <http://www.dhsspsni.gov.uk/niaic>  
Further information about **SABS** can be found at <http://sabs.dhsspsni.gov.uk/>

**Scotland**

Enquiries and adverse incident reports in Scotland should be addressed to:  
Incident Reporting and Investigation Centre  
NHS National Services Scotland, Gyle Square, 1 South Gyle Crescent, Edinburgh EH12 9EB  
Tel: 0131 275 7575 Fax: 0131 314 0722 Email: [nss.irc@nhs.net](mailto:nss.irc@nhs.net)  
<http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-irc/>

**Wales**

Enquiries in Wales should be addressed to:  
Improving Patient Safety Team  
Medical Directorate, Welsh Government, Cathays Park, Cardiff CF10 3NQ  
Tel: 029 2082 3922 Email: [Haz-Aic@wales.gsi.gov.uk](mailto:Haz-Aic@wales.gsi.gov.uk)