

June 25, 2013

## **URGENT PRODUCT CORRECTION NOTIFICATION** **Inspection of ORTHO BioVue® System Cassettes**

Dear Customer,

Ortho Clinical Diagnostics (OCD) has identified isolated occurrences of improperly positioned cassette labels for the products listed in the enclosure. Our records indicate that you were shipped some of the *potentially* affected products. The purpose of this notification is to inform you of the issue and provide instructions for the inspection of the potentially affected product remaining in your facility.

### **Background Information**

---

An investigation confirmed that the cause of the issue occurred on one of our three manufacturing lines following a particular sequence of events. The investigation concluded that the occurrence rate of the issue is very low.

The use of an affected multi-reagent cassette may lead to false negative or false positive results causing a potential misclassification of the patient or donor blood groups or incorrect antibody detection results. There is no risk associated with the use of an affected single-reagent cassette since all wells contain the same reagent.

### **Required Actions**

---

- Inspect all cassettes from the potentially affected lots remaining in your facility prior to use. Refer to the enclosed ORTHO BioVue® System Cassette Inspection Procedure for visual detailed instructions.
- Do not use cassettes with an incorrectly positioned label. Contact your Customer Technical Support representatives for to report the issue and request assistance. Discard the affected cassette in accordance with your local regulations.
- Complete and return the Confirmation of Receipt form no later than **July 2, 2013**.
- Post this notification and inspection procedure in your facility.
- Forward this notification if you have provided this product outside of your facility

### **Resolution**

---

We have identified the root cause of this issue and we have implemented corrective actions to mitigate a reoccurrence.

We apologize for the inconvenience this will cause your laboratory. We have anticipated some questions you may have in the following Question and Answers section. If you have any additional questions, please contact Customer Technical Services at 0800 895 963.

Sincerely,

Peter Clements  
Clinical Laboratories Marketing Manager Europe

Enclosures: Potentially Affected ORTHO BioVue® System Cassettes  
ORTHO BioVue® System Cassette Inspection Procedure

## Questions & Answers

**1. Is there any impact to previously reported results if I used a cassette from one of the affected lots?**

Our assessment determined that the likelihood an incorrect result obtained using an affected cassette is remote. This assessment is based on the following:

- Low probability of the occurrence of the issue
- Presence of a Control reagent in most blood grouping and phenotyping cassettes
- Use of Quality Control samples
- Laboratory procedures and/or specific regulations that require samples to be tested in duplicate or results are compared to previously obtained results

If you suspect that erroneous result occurred, please consult with your medical director and report the occurrence to our Customer Technical Service representatives.

**2. Are all lots affected by this issue?**

No, only the lots listed in the enclosure are subject to the inspection. Our investigation confirmed that the cause of the issue occurred on one of our three manufacturing lines following a particular sequence of events; not all types of cassettes are manufactured on this line. The investigation concluded that the occurrence rate of the issue is very low.

**3. How can I determine the proper label position if I have an affected lot?**

Follow the instructions provided in the enclosed ORTHO BioVue<sup>®</sup> System Cassette Inspection Procedure. The procedure sheet contains examples of both a properly and improperly positioned cassette label.

Performing the inspection prior to use will help to mitigate the occurrence and allow the use of the product.

**4. What action should I take if I identify an improperly positioned label cassette?**

If you identify an affected cassette in your facility, do not use the cassette. We advise that you contact your Customer Technical Support representatives to report the issue and request assistance. Discard the affected cassette in accordance with your local regulations.

**5. Does OCD recommend performing the inspection all at once?**

We recommend that you perform the inspection *prior to use* to minimize cassette handling and reduce the likelihood of other potential defects caused by excess handling.

**6. What is OCD doing to prevent this issue from reoccurring on other lots?**

We have identified the root cause of this issue and we have implemented preventative and corrective actions to mitigate a reoccurrence.

**Confirmation of Receipt - Important Response Required**

**URGENT PRODUCT CORRECTION NOTIFICATION**  
**Inspection of ORTHO BioVue® System Cassettes**

So that we can complete our records, please return this form to us no later than **July 2, 2013**.

**FAX TO:** *Shazia Younis*

**FAX:** *01494 658605*

**Section I: Confirmation**

I received the Urgent Product Correction Notification (Ref. CL13-198a) and understand that I must inspect all potentially affected cassettes (listed in the enclosure) for proper label orientation prior to use. If I identify an affected cassette, I am advised to not use it and contact my Customer Technical Support representatives to report the issue.

*Your signature provides confirmation that you have received and understood this notification.*

Your Name: \_\_\_\_\_

Job Title (optional): \_\_\_\_\_

Signed\*: \_\_\_\_\_

Date: \_\_\_\_\_

Fax Number: \_\_\_\_\_

Telephone Number: \_\_\_\_\_

J Number: \_\_\_\_\_

Institution: \_\_\_\_\_

*Your comments are always welcome:*

**Section II – Verification of your Name and Address**

Verify your name and mailing address:

Please complete this section if your name and/or mailing address have changed:

Institution / Contact Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State/Province: \_\_\_\_\_ Zip/Postal Code: \_\_\_\_\_

Telephone: \_\_\_\_\_ FAX: \_\_\_\_\_


# ORTHO BioVue<sup>®</sup> System Cassette Inspection Procedure

Step	Action
1.	Perform a visual inspection of the affected lots of ORTHO BioVue <sup>®</sup> System Cassettes prior to use. (Refer to CL13-198_Affected Products)
2.	Position the cassette so that the <b>Orientation Mark</b> is located on your right-hand side.
3.	<ul style="list-style-type: none"> <li>• If the <b>Name of the Reagent/Cassette</b> is facing you as shown in Figure 1a, it is acceptable to use the cassette for testing.</li> <li>• If the <b>Barcode</b> is facing you as shown in Figure 2b, do <b>not</b> use the cassette.</li> </ul>
4.	Do not use ORTHO BioVue <sup>®</sup> System Cassettes with an incorrectly positioned label (Figure 2b). Contact your Customer Technical Support representatives for to report the issue and request assistance at <b>insert appropriate number</b> .

## Properly Positioned Label

(Front View - Figure 1a)

(Back View -Figure 1b)



## Improperly Positioned Labels

(Front View - Figure 2a)

(Back View - Figure 2b)



**Affected Lots of ORTHO BioVue® System Cassettes Subject to Inspection**

<b>Product Code</b>	<b>Product Name</b>	<b>Lot No.</b>	<b>Expiry Date</b>
6904485	ABD Confirmation Cassette	NDC111A	7/4/2013
		NDC112A	10/6/2013
707135	ABD Confirmation Cassette	ACC515A	6/25/2013
		ACC516A	7/9/2013
		ACC517A	8/3/2013
		ACC518A	8/15/2013
		ACC519A	8/28/2013
		ACC520A	9/8/2013
		ACC521A	9/15/2013
		ACC522A	9/27/2013
		ACC523A	10/11/2013
		ACC524A	10/22/2013
		ACC526A	11/3/2013
		ACC527A	11/13/2013
		ACC528A	12/4/2013
		ACC529A	12/14/2013
		ACC530A	12/18/2013
		ACC531A	12/24/2013
		ACC531B	12/27/2013
		ACC532A	1/4/2014
ACC533A	1/9/2014		
ACC534A	1/11/2014		
ACC535A	1/16/2014		
ACC536A	1/30/2014		
707250	Rh/K Cassette	RHP332A	9/21/2013
		RHP335A	11/5/2013
		RHP337A	11/21/2013
		RHP338A	12/1/2013
		RHP339A	12/6/2013
		RHP343A	1/6/2014
		RHP345A	2/8/2014
707255	Rh-hr Cassette	BRH235A	8/4/2013
		BRH236A	9/10/2013
		BRH237A	9/28/2013

**Affected Lots of ORTHO BioVue® System Cassettes Subject to Inspection**

<b>Product Code</b>	<b>Product Name</b>	<b>Lot No.</b>	<b>Expiry Date</b>
<b>707280</b>	<b>Rh/K Cassette</b>	RHP327A	6/21/2013
		RHP330A	7/16/2013
		RHP332A	9/21/2013
		RHP335A	11/5/2013
		RHP337A	11/21/2013
		RHP338A	12/1/2013
		RHP339A	12/6/2013
		RHP340A	12/12/2013
<b>707119</b>	<b>ABO-DD Grouping Cassette</b>	ADD195A	7/17/2013
		ADD196A	8/27/2013
		ADD197A	10/15/2013
		ADD198A	11/5/2013
		ADD199A	12/29/2013
		ADD200A	1/15/2014
<b>707150</b>	<b>ABO-Rh Grouping Cassette</b>	ABE254A	7/14/2013
		ABE257B	8/25/2013
		ABE259A	10/8/2013
		ABE260A	10/16/2013
		ABE262A	11/15/2013
		ABE263A	12/11/2013
		ABE265A	12/21/2013
		ABE266A	12/24/2013
		ABE267A	1/7/2014
		ABE268A	1/10/2014
		ABE271A	2/18/2014
		ABE272A	2/22/2014

**Affected Lots of ORTHO BioVue® System Cassettes Subject to Inspection**

<b>Product Code</b>	<b>Product Name</b>	<b>Lot No.</b>	<b>Expiry Date</b>
<b>707190</b>	<b>ABO-Rh Grouping Cassette</b>	<b>ABE254A</b>	7/14/2013
		<b>ABE257A</b>	8/18/2013
		<b>ABE258A</b>	9/5/2013
		<b>ABE259A</b>	10/8/2013
		<b>ABE260A</b>	10/16/2013
		<b>ABE262A</b>	11/15/2013
		<b>ABE263A</b>	12/11/2013
		<b>ABE265A</b>	12/21/2013
		<b>ABE266A</b>	12/24/2013
		<b>ABE267A</b>	1/7/2014
		<b>ABE268A</b>	1/10/2014
		<b>ABE271A</b>	2/18/2014
		<b>ABE272A</b>	2/22/2014
<b>707310</b>	<b>Anti-Human Globulin Neutral Solution (Poly/Neutral Cassette)</b>	<b>PLN285A</b>	12/13/2013