

# Medical Device Alert

Ref: MDA/2013/013 Issued: 19 March 2013 at 14:00

## Device

Mitral Valve Repair System.

MitraClip Clip Delivery System.

Manufactured by Abbott  
Vascular.




Problem	Action
<p>Risk of death or serious harm to the patient.</p> <p>If the Actuator Knob on the Clip Delivery System (CDS) is turned in the wrong direction, this can prevent successful deployment of the clip, leading to the need for open surgical repair.</p>	<ul style="list-style-type: none"> <li>Follow the manufacturer's guidance given in the device's Instructions for Use (IFU) and in their recent <a href="#">Field Safety Notice (FSN)</a> with regards to the correct operation of the CDS Actuator Knob.</li> <li>Observe the blue arrow on the Actuator Knob, which is present on newer devices, indicating the direction in which it should be turned.</li> </ul>
Action by	
<p>Interventional Cardiologists Cardiothoracic Surgeons</p>	
CAS deadlines	Contact
<p>Action underway: 28 March 2013</p> <p>Action complete: 12 April 2013</p> <p><b>Note: These deadlines are for systems to be in place to implement actions</b></p>	<p>Anna Watson Abbott Vascular Tel: +44 (0) 1628 774977 Fax: +44 (0) 800 328 0868 Email: <a href="mailto:maivregulatory@av.abbott.com">maivregulatory@av.abbott.com</a></p>

## Problem


The MHRA has been notified by the manufacturer that there have been four reports (including one UK report), of the Actuator Knob of the Clip Delivery System (CDS) being turned in the wrong direction (i.e. clockwise instead of anti-clockwise) when attempting to disengage and deploy the clip from the CDS. Turning the Actuator Knob incorrectly could lead to a fracture of the CDS internal assembly, leading to an inability to control the release of the clip from the CDS. Since the clip may already have a firm grasp on the mitral valve at this point, open surgery may be required to disengage the clip from the CDS.

The UK incident led to the need for emergency open surgical repair. The bulk of the CDS could only be removed from the clip using surgical bolt cutters and the remaining shaft of the CDS unscrewed from the clip using forceps. Unfortunately, the patient died one week later, following renal complications.

In order to help prevent further incidents of this type, the manufacturer issued a FSN on 27 February 2013 providing advice. In addition, new devices will have a blue directional arrow on the Actuator Knob to help operators visualise the correct direction of turn. The manufacturer has also confirmed to the MHRA that they intend to make a number of design changes which should ultimately make it impossible to turn the Actuator Knob in the wrong direction. The implementation of the final stages of the design changes will not, however, be completed for up to two years.



From Figure 3 of the IFU,  
Identifying the Actuator Knob



## Distribution

This MDA has been sent to:

- NHS trusts in England (Chief Executives)
- Care Quality Commission (CQC) (Headquarters) for information
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Equipment Co-ordinators)
- NHS boards and trusts in Wales (Chief Executives)
- Primary care trusts in England (Chief Executives)

### **Onward distribution**

Please bring this notice to the attention of relevant employees in your establishment. Below is a suggested list of recipients.

#### **Trusts**

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- Cardiac surgical assistants
- Cardiologists
- Cardiothoracic surgeons
- Clinical governance leads
- Medical directors
- Nursing executive directors
- Supplies departments
- Theatre managers

### **Independent distribution**

#### **Establishments registered with the Care Quality Commission (CQC) (England only)**

This alert should be read by:

- Private hospitals

Please note: CQC and OFSTED do not distribute these alerts. Independent healthcare providers and social care providers can sign up to receive MDAs directly from the Department of Health's Central Alerting System (CAS) by sending an email to: [safetyalerts@dh.gsi.gov.uk](mailto:safetyalerts@dh.gsi.gov.uk) and requesting this facility.

## **Contacts**

#### **Manufacturer's Authorised Representative**

Jeremy Tinkler

MedPass International Ltd

Tel: +44 (0) 1452 619 222

Fax: +44 (0) 1452 619 222

Email: [Medpass.ar@medpass.org](mailto:Medpass.ar@medpass.org)

#### **Local contact for Manufacturer**

Anna Watson

Abbott Vascular

Tel: +44 (0) 1628 774977

Fax: +44 (0) 800 328 0868

Email: [maiavregulatory@av.abbott.com](mailto:maiavregulatory@av.abbott.com)

## **England**

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2013/013** or **2013/002/022/081/015**

#### **Technical aspects**

Alexander McLaren / Bayode Adisa

Medicines & Healthcare products Regulatory Agency

Floor 4

151 Buckingham Palace Road

London SW1W 9SZ

Tel: 020 3080 7292 / 7223

Fax: 020 8754 3965

Email: [alexander.mclaren@mhra.gsi.gov.uk](mailto:alexander.mclaren@mhra.gsi.gov.uk)

[bayode.adisa@mhra.gsi.gov.uk](mailto:bayode.adisa@mhra.gsi.gov.uk)

**Clinical aspects**

Nicola Lennard  
Medicines & Healthcare products Regulatory Agency  
Floor 4  
151 Buckingham Palace Road  
London SW1W 9SZ  
Tel: 020 3080 7126  
Fax: 020 8754 3965  
Email: [nicola.lennard@mhra.gsi.gov.uk](mailto:nicola.lennard@mhra.gsi.gov.uk)

**How to report adverse incidents**

Please report via our website <http://www.mhra.gov.uk>  
Further information about **CAS** can be found at <https://www.cas.dh.gov.uk/Home.aspx>

**Northern Ireland**

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.  
Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre  
Health Estates Investment Group  
Room 17  
Annex 6  
Castle Buildings  
Stormont Estate  
Dundonald BT4 3SQ  
Tel: 02890 523 704  
Fax: 02890 523 900  
Email: [NIAIC@dhsspsni.gov.uk](mailto:NIAIC@dhsspsni.gov.uk)  
<http://www.dhsspsni.gov.uk/index/hea/niaic.htm>

**How to report adverse incidents in Northern Ireland**

Please report directly to NIAIC, further information can be found on our website <http://www.dhsspsni.gov.uk/niaic>  
Further information about **SABS** can be found at <http://sabs.dhsspsni.gov.uk/>

**Scotland**

Enquiries and adverse incident reports in Scotland should be addressed to:  
Incident Reporting and Investigation Centre  
Health Facilities Scotland  
NHS National Services Scotland  
Gyle Square  
1 South Gyle Crescent  
Edinburgh EH12 9EB  
Tel: 0131 275 7575  
Fax: 0131 314 0722  
Email: [nss.irc@nhs.net](mailto:nss.irc@nhs.net)  
<http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-irc/>

## Wales

Enquiries in Wales should be addressed to:

Improving Patient Safety Team

Medical Directorate

Welsh Government

Cathays Park

Cardiff CF10 3NQ

Tel: 029 2082 3922

Email: [Haz-Aic@wales.gsi.gov.uk](mailto:Haz-Aic@wales.gsi.gov.uk)

MHRA is an executive agency of the Department of Health

© Crown Copyright 2013

Addressees may take copies for distribution within their own organisations