

# Medical Device Alert

Ref: MDA/2012/031 Issued: 06 June 2012 at 12:00

## Device

Anaesthetic machines and associated devices.

All manufacturers.

## Problem

The Association of Anaesthetists of Great Britain and Ireland (AAGBI) report that 'failure to check the anaesthetic machine and/or the breathing system features as a major contributory factor in many anaesthetic misadventures, including some that have resulted in hypoxic brain damage or death'.

## Action by

Anaesthetists, operating department practitioners, anaesthetic nurses and all staff using anaesthetic machines and associated devices.

## CAS deadlines

Action underway: 18 June 2012

Action complete: 25 June 2012

## Action

- Be aware that the AAGBI has published an updated version of '[Checking Anaesthetic Equipment](#)' endorsed by all UK Chief Medical Officers and the Royal College of Anaesthetists. The new 'Checklist for Anaesthetic Equipment 2012' is attached to this alert.
- Always follow the manufacturer's instructions for use.
- Check anaesthetic machines and all breathing system components before use.

## Distribution

This MDA has been sent to:

- NHS trusts in England (Chief Executives)
- Care Quality Commission (CQC) (Headquarters) for information
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Equipment Co-ordinators)
- Local authorities in Scotland (Equipment Co-ordinators)
- NHS boards and trusts in Wales (Chief Executives)

### Onward distribution

Please bring this notice to the attention of relevant employees in your establishment. Below is a suggested list of recipients.

### Trusts

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- A&E departments
- Anaesthetic nursing staff
- Anaesthetists
- Clinical governance leads
- Day surgery units
- Maternity units
- Medical directors
- Nursing executive directors
- Operating department practitioners
- Risk managers
- Theatre recovery suite
- Theatres

### Independent distribution

#### Establishments registered with the Care Quality Commission (CQC) (England only)

This alert should be read by:

- Hospitals in the independent sector
- Independent treatment centres

Please note: CQC and OFSTED do not distribute these alerts. Independent healthcare providers and social care providers can sign up to receive MDAs directly from the Department of Health's Central Alerting System (CAS) by sending an email to: [safetyalerts@dh.gsi.gov.uk](mailto:safetyalerts@dh.gsi.gov.uk) and requesting this facility.

## England

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2012/031** or **2012/005/022/291/017**

### Technical aspects

Douglas McIvor or Elke Kerwick  
Medicines & Healthcare products Regulatory Agency  
Floor 4  
151 Buckingham Palace Road  
London SW1W 9SZ

Tel: 020 3080 7193 / 6826

Fax: 020 8754 3965

Email: [douglas.mcivor@mhra.gsi.gov.uk](mailto:douglas.mcivor@mhra.gsi.gov.uk)  
[elke.kerwick@mhra.gsi.gov.uk](mailto:elke.kerwick@mhra.gsi.gov.uk)

**Clinical aspects**

Dr Nicola Lennard  
Medicines & Healthcare products Regulatory Agency  
Floor 4  
151 Buckingham Palace Road  
London SW1W 9SZ  
Tel: 020 3080 7126  
Fax: 020 8754 3965  
Email: [nicola.lennard@mhra.gsi.gov.uk](mailto:nicola.lennard@mhra.gsi.gov.uk)

**How to report adverse incidents**

Please report via our website <http://www.mhra.gov.uk>  
Further information about **CAS** can be found at <https://www.cas.dh.gov.uk/Home.aspx>

**Northern Ireland**

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.  
Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre  
Health Estates Investment Group  
Room 17  
Annex 6  
Castle Buildings  
Stormont Estate  
Dundonald BT4 3SQ  
Tel: 02890 523 704  
Fax: 02890 523 900  
Email: [NIAIC@dhsspsni.gov.uk](mailto:NIAIC@dhsspsni.gov.uk)  
<http://www.dhsspsni.gov.uk/index/hea/niaic.htm>

**How to report adverse incidents in Northern Ireland**

Please report directly to NIAIC, further information can be found on our website <http://www.dhsspsni.gov.uk/niaic>  
Further information about **SABS** can be found at <http://sabs.dhsspsni.gov.uk/>

**Scotland**

Enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre  
Health Facilities Scotland  
NHS National Services Scotland  
Gyle Square  
1 South Gyle Crescent  
Edinburgh EH12 9EB  
Tel: 0131 275 7575  
Fax: 0131 314 0722  
Email: [nss.irc@nhs.net](mailto:nss.irc@nhs.net)  
<http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-irc/>

## Wales

Enquiries in Wales should be addressed to:

Dr Chris Jones

Senior Medical Officer

Welsh Assembly Government

Cathays Park

Cardiff CF10 3NQ

Tel: 029 2082 3922

Email: [Haz-Aic@wales.gsi.gov.uk](mailto:Haz-Aic@wales.gsi.gov.uk)

MHRA is an executive agency of the Department of Health

© Crown Copyright 2012

Addressees may take copies for distribution within their own organisations

**Appendix**

**Checklist for Anaesthetic Equipment 2012**

**AAGBI Safety Guideline**



**Checks at the start of every operating session**  
**Do not use this equipment unless you have been trained**

**Check self-inflating bag available**

**Perform manufacturer's (automatic) machine check**

<b>Power supply</b>	<ul style="list-style-type: none"> <li>• Plugged in</li> <li>• Switched on</li> <li>• Back-up battery charged</li> </ul>
<b>Gas supplies and suction</b>	<ul style="list-style-type: none"> <li>• Gas and vacuum pipelines – ‘tug test’</li> <li>• Cylinders filled and turned off</li> <li>• Flowmeters working (if applicable)</li> <li>• Hypoxic guard working</li> <li>• Oxygen flush working</li> <li>• Suction clean and working</li> </ul>
<b>Breathing system</b>	<ul style="list-style-type: none"> <li>• Whole system patent and leak free using ‘two-bag’ test</li> <li>• Vaporisers – fitted correctly, filled, leak free, plugged in (if necessary)</li> <li>• Soda lime - colour checked</li> <li>• Alternative systems (Bain, T-piece) – checked</li> <li>• Correct gas outlet selected</li> </ul>
<b>Ventilator</b>	<ul style="list-style-type: none"> <li>• Working and configured correctly</li> </ul>
<b>Scavenging</b>	<ul style="list-style-type: none"> <li>• Working and configured correctly</li> </ul>
<b>Monitors</b>	<ul style="list-style-type: none"> <li>• Working and configured correctly</li> <li>• Alarms limits and volumes set</li> </ul>
<b>Airway equipment</b>	<ul style="list-style-type: none"> <li>• Full range required, working, with spares</li> </ul>

**RECORD THIS CHECK IN THE PATIENT RECORD**

<b>Don't Forget!</b>	<ul style="list-style-type: none"> <li>• Self-inflating bag</li> <li>• Common gas outlet</li> <li>• Difficult airway equipment</li> <li>• Resuscitation equipment</li> <li>• TIVA and/or other infusion equipment</li> </ul>
----------------------	--

This guideline is not a standard of medical care. The ultimate judgement with regard to a particular clinical procedure or treatment plan must be made by the clinician in the light of the clinical data presented and the diagnostic and treatment options available.  
 © The Association of Anaesthetists of Great Britain & Ireland 2012

## CHECKS BEFORE EACH CASE

### Breathing system

Whole system patent and leak free using 'two-bag' test  
 Vaporisers – fitted correctly, filled, leak free, plugged in (if necessary)  
 Alternative systems (Bain, T-piece) – checked  
 Correct gas outlet selected

### Ventilator

Working and configured correctly

### Airway equipment

Full range required, working, with spares

### Suction

Clean and working

## THE TWO-BAG TEST

**A two-bag test should be performed after the breathing system, vaporisers and ventilator have been checked individually**

- i. Attach the patient end of the breathing system (including angle piece and filter) to a test lung or bag.
- ii. Set the fresh gas flow to 5 l.min<sup>-1</sup> and ventilate manually. Check the whole breathing system is patent and the unidirectional valves are moving. Check the function of the APL valve by squeezing both bags.
- iii. Turn on the ventilator to ventilate the test lung. Turn off the fresh gas flow, or reduce to a minimum. Open and close each vaporiser in turn. There should be no loss of volume in the system.

This checklist is an abbreviated version of the publication by the Association of Anaesthetists of Great Britain and Ireland 'Checking Anaesthesia Equipment 2012'. It was originally published in *Anaesthesia*.  
 (Endorsed by the Chief Medical Officers)

If you wish to refer to this guideline, please use the following reference: Checklist for anaesthetic equipment 2012. *Anaesthesia* 2012; **66**: pages 662–63. <http://onlinelibrary.wiley.com/doi/10.1111/j.1365-2044.2012.07163.x/abstract>