

Medical Device Alert

Ref: MDA/2011/051 Issued: 19 May 2011 at 14:00

Device

Implantable cardioverter defibrillator (ICD) lead.

Sprint Fidelis, model numbers: 6930, 6931, 6948 and 6949

Manufactured by Medtronic.

Problem

Risk of inappropriate shocks and loss of therapy.

Increased risk that the high voltage conductor will fracture if the pace-sense conductor has previously fractured in the same lead.

Action by

All cardiologists and cardiac physiologists who manage patients implanted with ICD leads.

Action

Identify patients where:

- 1) the pace-sense conductor has fractured,
- 2) a new pace-sense lead has been inserted and
- 3) the high voltage conductor remains in use.

For these patients, consider the need to implant a new high voltage conductor at elective device change, taking into account the age and clinical circumstances of the patient.

If these patients are at high risk of sudden cardiac death in the event of high voltage conductor failure, then consider elective lead replacement before device change.

Where a new conductor fracture occurs, implant a new high voltage lead.

Where a lead is functioning normally, continue with existing patient management recommendations and use of Medtronic Lead Integrity Alert (LIA). This should be done taking into account the recent update in long term lead survival data presented in the [Field Safety Notice of April 2011](#).

CAS deadlines

Action underway: 26 May 2011

Action complete: 22 August 2011

Contact

Manufacturer

Mrs Lezlie Bridge
Medtronic

Tel: 01923 212 213

Fax: 01923 225 273

Email: lezlie.j.bridge@medtronic.com

Problem

In March 2009 Medtronic issued a [Field Safety Notice](#) (FSN) advising clinicians that one option for dealing with this problem was to implant a new pace-sense lead while maintaining the use of the high voltage conductor.

However, Medtronic estimates that the survival for the Sprint Fidelis high voltage conductor is 77.7% at 21 months if a pace-sense conductor fracture has previously occurred. Therefore, Medtronic no longer recommends this course of action and has issued another [FSN \(April 2011\)](#). All other patient management recommendations in the March 2009 FSN remain unchanged.

Distribution

This MDA has been sent to:

- NHS trusts in England (Chief Executives)
- Care Quality Commission (CQC) (Headquarters) for information
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Equipment Co-ordinators)
- NHS boards and trusts in Wales (Chief Executives)
- Primary care trusts in England (Chief Executives)

Onward distribution

Please bring this notice to the attention of relevant employees in your establishment. Below is a suggested list of recipients.

Trusts

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- A&E departments
- Arrhythmia nurses
- Cardiac laboratory technicians
- Cardiac pacing technicians
- Cardiologists
- Cardiology departments
- Cardiology nurses
- Cardiology, directors of
- Clinical governance leads
- Coronary care departments
- Coronary care nurses
- Medical directors
- Risk managers

Independent distribution

Establishments registered with the Care Quality Commission (CQC) (England only)

This alert should be read by:

- Independent treatment centres

Please note: CQC and OFSTED do not distribute these alerts. Independent healthcare providers and social care providers can sign up to receive MDAs directly from the Department of Health's Central Alerting System (CAS) by sending an email to: safetyalerts@dh.gsi.gov.uk and requesting this facility.

Contacts

Manufacturer

Mrs Lezlie Bridge
Medtronic Limited
Building 9
Croxley Green Business Park
Hatters Lane
Watford
Hertfordshire
WD18 8WW

Tel: 01923 212 213

Fax: 01923 225 273

Email: lezlie.j.bridge@medtronic.com

England

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2011/0051** or **2011/004/004/291/016**

Technical aspects

Ms Sam Baxter or Mr Guido Fumagalli
Medicines & Healthcare products Regulatory Agency
Floor 4
151 Buckingham Palace Road
London SW1W 9SZ

Tel: 020 3080 7242 / 7144

Fax: 020 8754 3965

Email: sam.baxter@mhra.gsi.gov.uk
guido.fumagalli@mhra.gsi.gov.uk

Clinical aspects

Dr Nicola Lennard
Medicines & Healthcare products Regulatory Agency
Floor 4
151 Buckingham Palace Road
London SW1W 9SZ

Tel: 020 3080 7126

Fax: 020 8754 3965

Email: nicola.lennard@mhra.gsi.gov.uk

How to report adverse incidents

Please report via our website <http://www.mhra.gov.uk>

Further information about **CAS** can be found at <https://www.cas.dh.gov.uk/Home.aspx>

Northern Ireland

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.

Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre

Health Estates Investment Group

Room 17

Annex 6

Castle Buildings

Stormont Estate

Dundonald BT4 3SQ

Tel: 02890 523 704

Fax: 02890 523 900

Email: NIAIC@dhsspsni.gov.uk

<http://www.dhsspsni.gov.uk/index/hea/niaic.htm>

How to report adverse incidents in Northern Ireland

Please report directly to NIAIC, further information can be found on our website <http://www.dhsspsni.gov.uk/niaic>

Further information about **SABS** can be found at <http://sabs.dhsspsni.gov.uk/>

Scotland

Enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre

Health Facilities Scotland

NHS National Services Scotland

Gyle Square

1 South Gyle Crescent

Edinburgh EH12 9EB

Tel: 0131 275 7575

Fax: 0131 314 0722

Email: nss.irc@nhs.net

<http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-irc/>

Wales

Enquiries in Wales should be addressed to:

Dr Sara Hayes

Senior Medical Officer

Medical Device Alerts

Welsh Assembly Government

Cathays Park

Cardiff CF10 3NQ

Tel: 029 2082 3922

Email: Haz-Aic@wales.gsi.gov.uk

MHRA is an executive agency of the Department of Health

© Crown Copyright 2011

Addressees may take copies for distribution within their own organisations