Medical Device Alert



Action

Ref: MDA/2010/073 Issued: 20 September 2010 at 14:30

Device

Intravenous (IV) extension sets

with multiple ports: all brands.



Problem	Action
Risk of back-tracking when an IV line has multiple access ports.	 Consider if an IV extension set with multiple ports is appropriate for your intended use and be aware of alternatives. Where appropriate, consider using IV lines with one-way valves to prevent back-tracking (examples of one-way valves are: check, non-return or anti-reflux valves or anti-siphon/anti free-flow valves).
This can lead to under-infusion or bolus delivery of IV drugs.	
Action by	
All medical and nursing staff.	
CAS deadlines	 Apply clamps (where available) to lines not in use.
Action underway: 04 October 2010	Be aware that needle-free connectors are not and will allow back tracking
Action complete: 01 November 2010	one-way valves and will allow back-tracking when connected to IV devices.
	This MDA supersedes MDA/2007/089.

Device

All IV lines with multiple ports that will allow fluid to flow in both directions:

- multi-lumen IV extension sets
- IV administration sets with side port/Y-site
- IV multi-way connectors such as stopcocks
- Y-connectors
- manifolds.

This issue is not limited to any particular manufacturer or model.

Problem

The MHRA continues to receive reports of incidents involving intravenous (IV) extension sets with multiple ports. These incidents have led to serious consequences e.g. patients receiving an inadvertent bolus of anaesthetic agent due to back-tracking, resulting in respiratory arrest.

When more than one IV line is connected to a single access point, back-tracking of fluids can occur if one line has no flow or a slower flow of fluid running through it. The fluid will back-track and take the path of least resistance. This can also occur during occlusions. If a downstream occlusion occurs when an infusion pump is being used, the occlusion alarm may not be activated.

Such back-tracking can lead to bolus delivery when: flow in the line is increased; or the line is flushed; or the occlusion is released.

Distribution

This MDA has been distributed to:

- NHS trusts in England (Chief Executives)
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Chief Executives)
- NHS boards and trusts in Wales (Chief Executives)
- Primary care trusts in England (Chief Executives)

Onward distribution

Please bring this notice to the attention of all who need to know or be aware of it. This may include distribution by:

Trusts to:

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- Accident and emergency departments
- All clinical departments
- All wards
- Ambulance services directors
- Ambulance staff
- Anaesthesia, directors of
- Chief pharmacists
- Clinical governance leads
- Hospital at home units
- Intensive care unitsIV nurse specialists
- Medical directors
- Neonatal nurse specialists
- Nursing executive directors
- Nutritional nurse specialists
- Outpatient clinics
- Outpatient theatre managers
- Palliative care teams
- Paramedics
- PharmacistsRadiographer superintendents
- Resuscitation officers and trainers
- Risk managers
- Supplies managers
- Theatre managers

Establishments registered with the Care Quality Commission (CQC) (England only) to:

The MHRA considers this information to be important to:

- Adult placement
- Care homes providing nursing care (adults)
- Hospices
- Hospitals in the independent sector
- Independent treatment centres
- Private medical practitioners

Primary care trusts to:

CAS liaison officers for onward distribution to all relevant staff including:

- Community children's nurses
- Community hospitals
- Community midwives
- Community nurses
- District nurses
- Nutritional nurse specialists
- Palliative care teams
- NHS walk in centres

England

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2010/073** or **2010/002/023/601/003**.

Technical aspects

Dr Louise Mulroy or Mrs Nicole Small Medicines & Healthcare products Regulatory Agency Market Towers 1 Nine Elms Lane London SW8 5NQ

Tel: 020 7084 3344 Fax: 020 7084 3565

Email: louise.mulroy@mhra.gsi.gov.uk nicole.small@mhra.gsi.gov.uk

Clinical aspects

Dr Susanne Ludgate Medicines & Healthcare products Regulatory Agency Market Towers 1 Nine Elms Lane London SW8 5NQ

Tel: 020 7084 3123 Fax: 020 7084 3111

Email: susanne.ludgate@mhra.gsi.gov.uk

How to report adverse incidents

Please report via our website http://www.mhra.gov.uk Further information about **CAS** can be found at https://www.cas.dh.gov.uk/Home.aspx

Northern Ireland

Alerts in Northern Ireland will continue to be distributed via the NI SABS system. Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre Health Estates Investment Group Room 17 Annex 6 Castle Buildings Stormont Estate Dundonald BT4 3SQ Tel: 02890 523 704

Fax: 02890 523 900

Email: NIAIC@dhsspsni.gov.uk http://www.dhsspsni.gov.uk/index/hea/niaic.htm

How to report adverse incidents in Northern Ireland

Please report directly to NIAIC, further information can be found on our website http://www.dhsspsni.gov.uk/niaic Further information about **SABS** can be found at http://sabs.dhsspsni.gov.uk/

Scotland

Enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre Health Facilities Scotland NHS National Services Scotland Gyle Square 1 South Gyle Crescent Edinburgh EH12 9EB

Tel: 0131 275 7575 Fax: 0131 314 0722

Email: nss.iric@nhs.net

http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-iric/

Wales

Enquiries in Wales should be addressed to:

Dr Sara Hayes Senior Medical Officer Medical Device Alerts Welsh Assembly Government Cathays Park Cardiff CF10 3NQ

Tel: 029 2082 3922

Email: Haz-Aic@wales.gsi.gov.uk

MHRA is an executive agency of the Department of Health © Crown Copyright 2010 Addressees may take copies for distribution within their own organisations