

Medical Device Alert

Action

Ref: MDA/2009/080 Issued: 09 December 2009 at 15:30

Device

Ultrasound transducer probes with an internal lumen used for taking transrectal prostate biopsies.

All manufacturers.

Problem	Action
<p>Potential onward transmission of abnormal prion protein, the putative infective agent in variant Creutzfeldt Jakob disease (vCJD), when ultrasound transducer probes with an internal lumen are used for taking transrectal prostate biopsies on men at risk of vCJD.</p> <p>The needle that passes through the internal lumen of a reusable ultrasound probe could transfer the abnormal prion protein from rectal lymphoid tissue, via the outside surface of the biopsy needle, onto the internal lumen of the ultrasound probe. Any attempt at subsequent decontamination of the internal lumen would not be successful due to the inability of any current decontamination process to reliably eliminate or destroy abnormal prion protein.</p>	<p>Review the advice given in the Advisory Committee on Dangerous Pathogens (ACDP) TSE Working Group's alert to urological surgeons 'Transrectal prostatic biopsy in men at risk of variant CJD', published on the Department of Health's website.</p> <p>The ACDP TSE Working Group and CJD Incidents Panel advise the following: 'Patients at risk of vCJD requiring transrectal prostatic biopsy should have the procedure performed by means of single use equipment that runs alongside (rather than through) the ultrasound probe. Where a unit does not have such equipment and intends to carry out a biopsy procedure on a patient at risk of vCJD, their options are as follows:</p> <ul style="list-style-type: none"> • To refer the patient to a unit offering an alternative technique that does not pose a risk of contaminating the internal channels with traces of biopsy tissue • To borrow the alternative equipment from another unit • To undertake the procedure with equipment that has internal biopsy channels, and then quarantine the reusable components of that equipment after decontamination. It must be accepted that this equipment would be unlikely to return to general use, except for dedicated re-use in the same patient.'
<h3>Action by</h3>	
<p>All staff who carry out prostate biopsies.</p>	
<h3>CAS deadlines</h3>	
<p>Action underway: 31 December 2009 Action complete: 26 February 2010</p>	

Device

All ultrasound transducer probes which require the passing of a biopsy needle through the internal lumen.

Action

From the ACDP TSE Working Group's alert to urological surgeons:

'The following patient groups have been notified of their increased risk of subclinical vCJD infection:

- people who have received blood from someone who went on to develop vCJD
- people who have given blood to someone who went on to develop vCJD
- people who have received blood from someone who has also given blood to a patient who went on to develop vCJD
- people who have had surgery using instruments that had been used on someone who developed vCJD
- people who have had a neurosurgical procedure, or an operation for a tumour or cyst of the spine, before August 1992
- people who have received an organ or tissue from a donor infected with vCJD or at increased risk of vCJD
- people who have been treated with certain UK sourced plasma products between 1980 and 2001.

It is important to note that new patient groups may be notified in the future of their increased risk of vCJD.'

The advice in the ACDP's alert does not apply in procedures on men who are not at risk of vCJD nor to the use of probes that are covered by latex in their entirety while guided prostatic biopsy is being performed via a separate single use external biopsy sheath.

Distribution

This MDA has been distributed to:

- NHS trusts in England (Chief Executives)
- Care Quality Commission (CQC) (Headquarters)
- Health Protection Agency (HPA) (Directors)
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Chief Executives)
- NHS boards and trusts in Wales (Chief Executives)
- Primary care trusts in England (Chief Executives)

Onward distribution

Please bring this notice to the attention of all who need to know or be aware of it. This may include distribution by:

Trusts to:

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- Clinical governance leads
- Day surgery units
- Decontamination leads
- Directors of infection prevention and control
- Directors of nursing
- General surgeons
- General surgical units, directors of
- Health and safety managers
- Infection control departments
- Infection control nurses
- Medical directors
- Medical physics departments
- Microbiologists
- Oncologists

- Outpatient clinics
- Outpatient theatre managers
- Outpatient theatre nurses
- Purchasing managers
- Radiographer superintendents
- Radiographers
- Radiologists
- Radiology departments
- Radiology directors
- Risk managers
- Sterile services departments
- Supplies managers
- Theatre managers
- Theatre nurses
- Theatres
- Urological surgeons
- Urological surgery, directors of
- Urology departments

Care Quality Commission (CQC) (England only) to:

The MHRA considers this information to be important to:

- Adult placement
- Clinics
- Hospitals in the independent sector
- Independent treatment centres
- Private medical practitioners

Health Protection Agency to:

Directors for onward distribution to:

- Consultants in communicable disease control
- Health protection nurses
- Regional epidemiologists
- Regional microbiologists

Primary care trusts to:

CAS liaison officers for onward distribution to all relevant staff including:

- Community hospitals
- Directors of infection prevention and control

England

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2009/080** or **2009/010/030/291/001**

TSEs and decontamination issues

Allan Hilderley
Medicines & Healthcare products Regulatory
Agency
Market Towers
1 Nine Elms Lane
London SW8 5NQ
Tel: 020 7084 3172
Fax: 020 7084 3106
E-mail: allan.hilderley@mhra.gsi.gov.uk

Technical issues

Richard Glover
Medicines & Healthcare products Regulatory
Agency
Market Towers
1 Nine Elms Lane
London SW8 5NQ
Tel: 020 7084 3245
Fax: 020 7084 3209
E-mail: richard.glover@mhra.gsi.gov.uk

Clinical aspects

Christopher Brittain
Medicines & Healthcare products Regulatory Agency
Market Towers
1 Nine Elms Lane
London SW8 5NQ
Tel: 020 7084 3126
Fax: 020 7084 3111
E-mail: christopher.brittain@mhra.gsi.gov.uk

How to report adverse incidents

Please report via our website <http://www.mhra.gov.uk>

Further information about **CAS** can be found at <https://www.cas.dh.gov.uk/Home.aspx>

Northern Ireland

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.

Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre (NIAIC)

Health Estates

Estate Policy Directorate

Stoney Road

Dundonald

Belfast

BT16 1US

Tel: 02890 523 704

Fax: 02890 523 900

E-mail: NIAIC@dhsspsni.gov.uk

<http://www.dhsspsni.gov.uk/index/hea/niaic.htm>

How to report adverse incidents in Northern Ireland

Please report directly to NIAIC, further information can be found on our website <http://www.dhsspsni.gov.uk/niaic>

Further information about **SABS** can be found at <http://sabs.dhsspsni.gov.uk/>

Scotland

Health Protection Scotland (HPS) advises that the advice given in the Advisory Committee on Dangerous Pathogens (ACDP) TSE Working Group's Alert titled 'Transrectal prostatic biopsy in men at risk of variant CJD' is applicable in Scotland.

Enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre

Health Facilities Scotland

NHS National Services Scotland

Gyle Square

1 South Gyle Crescent

Edinburgh

EH12 9EB

Tel: 0131 275 7575

Fax: 0131 314 0722

E-mail: nss.irc@nhs.net

<http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-irc/>

Wales

Dr Jane Ludlow
Senior Medical Officer
Medical Device Alerts
Welsh Assembly Government
Cathays Park
Cardiff
CF10 3NQ
Tel: 029 2082 3505 / 3922
E-mail: Haz-Aic@wales.gsi.gov.uk

MHRA is an executive agency of the Department of Health
© Crown Copyright 2009
Addressees may take copies for distribution within their own organisations