

Anti-Malarial Drug Testing: Tests, Reference Standards and Assays

SensaPharm Ltd, Sunderland UK

Enquiries to: neil@sensapharm.eu

+44 191 516 6969

Part funded by;



Medicines for Malaria Venture



Topics

- Reference Standards
- Validated Assays for Artemisinin
- Laboratory Tests
- In-Field Tests
- Metabolic profiling, plant characterisation and traceability for A₂S₂

Reference Standards

- Current supply of anti-malarial standards is limited and sparse
- Comprehensive portfolio of standards, fully characterised
- Contains 40 standards and 20 plant extracts
 - artemisinin and artemisinin derived pharmaceuticals
 - other anti-malarial pharmaceuticals
 - artemisinin derivatives
 - Useful plant derived metabolites
 - *A. annua* variety Artemis extracts and leaf samples

Reference Standards - Examples

ART & ART derived pharmaceuticals

- Artemisinin
- Arteether
- Artemether
- Artesunate

Other anti-malarial pharmaceuticals

- Lumefantrine
- Piperaquine phosphate
- Sulphamethoxypyrazine
- Trimethoprim

A. *annua* variety Artemis extracts

- Hexane- EtOAc
- R_{134A}
- Acetonitrile
- Toluene

Artemisinin derivatives

- Dihydroartemisinin
- Deoxyartemisinin

Useful plant derived metabolites

- Quinine, Quercetin, Coumarin
- Artemisinic acid, DHAA
- Arteannuin B, Artemisitene

Artemisia leaf samples

- *Artemisia abrotanum*
- *Artemisia absinthium*
- *Artemisia vulgaris*

Validated Assay for Artemisinin using LCMSMS



The need for an accredited assay for the detection and measurement of ART

- Current methodologies are not best suited to the detection and measurement of Artemisinin
- Current methods for ‘pure’ Artemisinin analysis could be improved.
- We need an assay that is suited for quantitative analysis of both ‘pure’ Artemisinin and Artemisinin present in leaf extracts with options for its derivatives and impurities.

Validated Assay


- We have now launched a validated LCMSMS assay service for artemisinin and related compounds
- BS 17025 accreditation assessment in November and is based upon ICH Q2A guidelines
- Available for inter laboratory calibration...not intended to replace existing assays, but to help provide consistency of analysis in the industry
- ...send us 1 in every 10 samples for calibration

Antibody based detection systems for Artemisinin


ELISA and Lateral Flow Devices




Overview



ELISA for Artemisinin;




- Limits of detection = 1mg/ml-50ng/ml
- Limits of quantification = 2-64µg/ml



- ELISA give excellent correlation with other analytical techniques
- Supplement to TLC, more reproducible
- Use in conjunction with validated LCMSMS

In Field Quantitation Kit

- For use with Lateral Flow Device and Device reader
- Configured for field use with all equipment an instructions included
- Looking for beta testing partners



SensaPharm

Principle Component Analysis
using Multivariate Data Analysis
(PCA) for metabolic profiling in
A. annua

The advantage of PCA:

- A mathematical technique based upon real data which allows identification of molecules responsible for key characteristics
- It is the only way to clearly identify the molecules (or groups of molecules) which are likely to relate to industry issues such as; agronomy (region of growth, pesticides, fertilizers, irrigation, sunlight and daylight effects, high yielding plants, extraction and crystallization efficiency, API impurities*, API stability, future processing
- This technique will give us the tools required for rational design of industry protocols.
- * Important implications for supply chain traceability from plant to API to final product

Data for PCA:

VARIABLES

O
B
S
E
R
V
A
T
I
O
N
S

Primary Obs	ClassID	Maturity	2.564"	DA2.659"	DA2.892"	DA3.048"	DA3.219"	DA3.40"	DAE3.50"	DAE3.653"	DA3.842"	DA4.0"	DAD14.1"	DAD14.303"	DA4.475"	DA4.739"	DA4.8"	DAD
S082-37	Iboaka 1045 Mat-4,2	4,2	26.86961	44.9733	52.41531	54.34677	51.66421	123.5335	0	58.88955	41.57663	0	0	631.0388	0	174.1241	741.97	
S082-38	Iboaka 1046 Mat-4,4	4,4	25.25257	43.72847	51.26069	0	142.8474	155.5588	0	0	246.9885	557.2141	0	0	76.8434	0	846.5731	
S082-39	Ambalaseva 3194 Mat-	Ambalaseva	26.86961	50.38239	82.44022	0	218.0575	100.0604	0	0	878.6672	0	0	78.72085	89.85306	0	805.6885	
S082-41	Antisrabe 1225 Mat-4	4	1093.804	0	611.4351	0	0	59.09922	518.197	335.384	0	0	78.61356	65.32234	80.53956	0	646.5154	
S082-42	Antisrabe 1226 Mat-4,9	4,9	1060.033	0	479.2135	0	111.7254	0	965.2594	0	0	88.69537	81.01097	0	226.5649	0	902.8701	
S082-43	Antisrabe 1021 Mat-3,6	3,6	953.4064	0	308.4754	195.9851	0	0	571.5104	500.549	0	0	50.18599	0	102.5747	0	518.682	
S082-44	Antisrabe 1023 Mat-3,1	3,1	857.1718	0	294.5936	176.8587	0	0	428.0541	536.7195	0	0	43.52091	0	78.47832	0	416.9307	
S082-45	Iboaka 1238 Mat-5,1	5,1	123.6765	69.77823	334.3663	0	136.9739	0	1093.075	0	0	70.13799	59.47911	0	176.0726	0	942.0013	
S082-46	Iboaka 1201 Mat-4,7	4,7	145.984	68.87493	379.0592	0	135.1142	0	842.2739	0	0	46.51926	58.01641	63.22187	75.76952	0	741.4986	
S082-47	Iboaka 1204 Mat-4,9	4,9	0	0	124.6906	161.2865	76.73167	0	316.0786	613.0707	0	0	0	129.1609	93.39915	0	731.0713	
Med ACN1	Mediplant	Med ACN1	177.3451	180.6178	105.9363	320.3188	86.90177	601.4019	0	199.4794	143.51	0	0	144.8257	0	0	0	
Med ACN2	Mediplant	Med ACN2	168.3885	181.7006	105.922	332.5286	94.85378	602.0241	0	282.0712	0	0	53.20174	159.8341	0	0	0	
Med ACN3	Mediplant	Med ACN3	163.7602	294.3558	0	338.8649	98.37358	597.8733	0	284.6949	50.27254	0	0	163.7376	0	0	0	


- PCA breaks down a large table of data into two smaller ones
- Plots of scores and loadings turn data into pictures
- Correlations among observations and variables are easily seen
- Information resides in *the correlation structure* of the data
- Interpretation of PCA scores shows patterns and trends



Data acquisition for PCA:

- We initially chose to use HPLC at 216 and 280 nm to establish the model
- Q-TOF MS under way...to identify **all** components in a plant and to analyze against relevant variables
- Studies in progress with FTIR as rapid analysis technique

Example analysis...this is for a particular localized region of growth, but we have also examined the compounds extracted by different solvents, from internationally sourced leaf samples and from different cultivars.



Yellow- Mediplant std.

Others- Same variety grown in different areas...they can be distinguished by molecular analysis

Example of PCA model Observations

- In proof of principle studies we have shown that the components of *A. annua* leaf from different geographical origins, and extracted using different solvents can be separated into distinct groups
- The differences in metabolic profiles which are impossible to interpret with the naked eye, can be identified using multivariate analysis e.g. to identify compounds related to crystallisation and yield
- The models developed may also be used for monitoring the batch-to-batch consistency of raw materials and for traceability in the supply chain (A2S2)

Summary

- **Our aim is to promote market normalisation and stability by providing a helpful portfolio of relevant products and services. These include;**
- A new portfolio of anti-malarial reference standards
- New in-field and laboratory based tests
- Validated assay services for artemisinin, several other related compounds and other anti-malarials
- Metabolic profiling and supply chain traceability using multivariate data analysis
- **Please contact us at: neil@sensapharm.eu**

