



Medical Device Alert

MDA/2019/033

Issued: 08 October 2019 at 14:00

Anaesthetic face masks – Specific Intersurgical Economy 22F taper connection may be oversized and leak or disconnect from the breathing circuit

Summary

Manufactured by Intersurgical – leaks or disconnection can result in insufficient oxygenation requiring medical intervention to avoid severe injury.

Action

- Check all stock for affected devices listed in the manufacturer's [Field Safety Notice \(FSN\)](#).
- Before use, check the connection between the taper and the mating device as described in the [FSN](#) and instructions for use.
- If you identify any problems with these devices, contact the manufacturer to arrange return and replacement of the devices.
- Complete the response form supplied in the [FSN](#) and return to Intersurgical, even if you no longer have affected devices in stock.
- Report suspected or actual adverse events involving these devices through your local incident reporting system and/or your national incident reporting authority as appropriate: [England](#), [Scotland](#), [Northern Ireland](#), [Wales](#). You should also report directly to manufacturers if your local or national systems do not.

Action by

All healthcare professionals who are responsible for or who use these devices.

Deadlines for actions

Actions underway: 15 October 2019

Actions complete: 29 October 2019

Medical Device Safety Officers (in England): ask the manufacturer to add you to their distribution list for field safety notices (FSNs). This is to help with reconciliation.

Remember: if your organisation receives an [FSN](#) from a manufacturer, always act on it. **Do not wait** for a communication from MHRA.

Manufacturer contacts

Intersurgical

Tel: 0118 9656362

Email: priority@intersurgical.co.uk

Distribution

If you are responsible for cascading these alerts in your organisation, these are our suggested distribution lists.

Trusts (NHS boards in Scotland)

CAS and NICAS liaison officers for onward distribution to all relevant staff including:

- A&E consultants
- A&E departments
- A&E directors
- A&E nurses
- Adult intensive care units
- Ambulance services directors
- Ambulance staff
- Anaesthesia, directors of
- Anaesthetic medical staff
- Anaesthetic nursing staff
- Anaesthetists
- Clinical perfusionists
- Day surgery units
- EBME departments
- Equipment stores
- Intensive care medical staff/paediatrics
- Intensive care nursing staff (adult)
- Intensive care nursing staff (paediatric)
- Intensive care units
- Intensive care, directors of
- Operating department practitioners
- Oral surgeons
- Outpatient theatre managers
- Outpatient theatre nurses
- Paediatric intensive care units
- Paediatric wards
- Paramedics
- Resuscitation officers and trainers
- Risk managers
- Special care baby units
- Supplies managers
- Theatre managers
- Theatre nurses
- Theatres

Independent distribution

Establishments registered with the Care Quality Commission (CQC) (England only)

- Hospitals in the independent sector
- Independent treatment centres
- Private medical practitioners

Please note: CQC and OFSTED do not distribute these alerts. Independent healthcare providers and social care providers can sign up to receive MDAs directly from the Central Alerting System (CAS) by sending an email to: safetyalerts@mhra.gov.uk and requesting this facility.

Enquiries

England

Send enquiries about this notice to MHRA, quoting reference number MDA/2019/033 or 2019/007/010/228/003.

Technical aspects

Ben Satchell or Enitan Taiwo, MHRA

Tel: 020 3080 6000

Email: DSS-TM@mhra.gov.uk

Clinical aspects

Devices Clinical Team, MHRA

Tel: 020 3080 7274

Email: dct@mhra.gov.uk

To report an adverse incident involving a medical device in England use the [Yellow Card reporting page](#)

Northern Ireland

Northern Ireland Adverse Incident Centre (NIAIC), CMO Group, Department of Health (Northern Ireland)

Tel: 028 9052 3868

Email: niaic@health-ni.gov.uk

To report an adverse incident involving a medical device in Northern Ireland use the [forms on the website](#).

Alerts in Northern Ireland are distributed via the [NICAS system](#).

Scotland

Incident Reporting and Investigation Centre (IRIC), Health Facilities Scotland, NHS National Services Scotland

Tel: 0131 275 7575

Email: nss.irc@nhs.net

To report an adverse incident involving a medical device in Scotland, [email IRIC](#) to request a webform account.

For more information, or if you can't access the webform, visit the website: [how to report an adverse incident](#)

Wales

Population Healthcare Division, Welsh Government

Tel: 03000 255278 / 03000 255510

Email: haz-aic@gov.wales

To report an adverse incident involving a medical device in Wales, use the [Yellow Card reporting page](#) and follow specific advice for reporting in Wales in [MDA/2004/054 \(Wales\)](#).

MHRA is a centre of the Medicines and Healthcare products Regulatory Agency, an executive agency of the Department of Health and Social Care.

© Crown Copyright 2019

Addressees may take copies for distribution within their own organisations