



HONEY: NMR Authenticity testing

Dr. Eric Jamin
Authenticity Competence Centre Manager
Eurofins Analytics France (Nantes)

Context



- Honey is one of the « top 10 » products regarding the risk of food fraud, according to an EU report dated 8/10/2013:
Draft report on the food crisis, fraud in the food chain and the control thereof 2013/2091 (INI)
- Fewer bee hive survival & lower yields worldwide, but yet global production increase, with large differences in price
- **The main fraud encountered in honey is sugar addition**
- **Other frauds may concern mislabelling of the botanical and/or geographic origin**
- This slide show aims to describe the state-of-the-art and Eurofins analytical offer, including recent developments

Our vision of honey market needs



 eurofins

Analyses must help to enforce regulations & protect brands but:



- Need for a high number of parameters
- Often limited time, and
- Always a limited budget



- Targeted analysis efficiency is limited, as fraudsters are always pushing limits



**Need to select the most efficient tests
& Include untargeted screening**

3

Optimal offer covering the main potential fraud risks



 eurofins

- **NMR profiling**
 - broad authenticity and integrity screening
- **Stable isotope testing via EA-IRMS**
 - detecting C4 sugars with an optimum detection limit
- **Pollen/sensorial/conductivity analysis**
 - confirming botanical/geographical origin and quality
- **LC-HRMS**
 - detecting some sugars with an optimum detection limit

NEW

4

Honey Expert laboratories



Our expert laboratories are ISO 17025 accredited*
for all above mentioned honey tests



Eurofins Authenticity competence centre (Nantes, France)

* accredited by COFRAC under registration n° 1-0287, details on www.cofrac.fr



Eurofins Food Integrity Control Service GmbH (Bremen, Germany)

* accredited by DAkkS under registration n° D-PL-20965-01-00, details on www.dakks.de

14/11/2019

Authenticity testing

5

NMR profiling High Resolution ^1H NMR



Observation domain:
from ppm to %




Key Strengths:

Broad observation of the whole matrix

High Reproducibility

Discrimination power



Eurofins initial publication

Food Chemistry 189 (2015) 60–66



Contents lists available at ScienceDirect

Food Chemistry

journal homepage: www.elsevier.com/locate/foodchem



— Lab 1
— Lab 2

Fast and global authenticity screening of honey using ^1H -NMR profiling



Marc Spiteri^a, Eric Jamin^{a,b}, Freddy Thomas^a, Agathe Rebours^a, Michèle Lees^a, Karyne M. Rogers^b, Douglas N. Rutledge^c

^a Eurofins Analytics France, Rue Pierre Adolphe Béhéret, B.P. 42301, F-44323 Nantes Cedex 3, France

^b National Isotope Centre, GNS Science, 30 Grandfield Road, Seaview, P.O. Box 31 312, Lower Hunt 5040, New Zealand

^c AgrobiosTech, UMR1145 Ingénierie Produits Alimentés, 16, rue Claude Bernard, 75005 Paris, France

Since 2018:
SOP aligned with Bruker
Honey consortium

6


NMR profiling Benefits of a non-targeted method



eurofins

- The $^1\text{H-NMR}$ spectrum is obtained on the whole matrix

Unique spectral fingerprint



- Detects markers and spectral effects of adulteration
- Allows to quantify key-quality parameters
- Allows to confirm botanical origins of mono-floral honey

Authenticity proof

Alternative to
conventional methods

7


NMR profiling Building a worldwide authentic honey database



eurofins

Reference samples : > 10000
>100 different botanical « families »

Collected from local producers
In more than 65 countries worldwide
over more than 15 years




The inclusion of samples in the statistical models requires:

- Availability of metadata (botanical & geographical origin)
- Compliance of other methods results
- Coherence with other samples of the same kind (95% confidence intervals)

8

**Analytical Method lifecycle:
a continuous improvement process
within flexible accreditation**



9

EDR: Auteur: Date:

**PRO-PTS: Proficiency Testing
for NMR PROFILING methods**



METHODS CONCERNED

- Wine NMR profiling
- Juice NMR profiling
- Honey NMR profiling
- Coffee NMR profiling
- Spices NMR profiling
- Oils NMR profiling

Since 2017: 3 times a year
Ca 20 participants (10 for honey)
Each lab is free to use its internal SOP


STATISTICAL ANALYSIS

- Calculation of z-scores using robust statistics
- According to the ISO/IUPAC/AOAC International Harmonised Protocol for Proficiency Testing of analytical laboratories.

- Includes non targeted analysis / data interpretation
- Cheks global agreement between laboratories

Stat	2019 Round 2	Product	Parameter
Min	0.0	Honey	Moisture
Max	2.2		
Target Mean	1.0		
Standard Deviation	0.5		
Reference SD	0.5		
Target SD	0.2		
Min	2.2		
Max	2.2		

Reference Code	Sample	Element
E1	2.0	-0.76
E7	2.0	-0.73
E1	2.0	-0.20
E8	2.0	-0.20
E2	2.0	-0.20
E3	2.0	-0.17
E4	2.0	-0.17
E5	2.0	-0.15
E6	2.0	-0.15
E9	2.0	-0.15
E10	2.0	-0.15
E11	2.0	-0.15
E12	2.0	-0.15



2019 Round 2		Product: Honey (moStk)			Sample code	19-28
Lab.n°	Address	Honeytype	France	Deviations from model		
E1	SC	Yes	Yes	/	/	/
E2	SC	Yes	Yes	/	/	/
E3	SC	Yes	Yes	/	/	/
E4	SC	Yes	Yes	Reference construction elevated	/	/
E5	SC	Yes	Yes	/	/	/
E6	SC	Yes	Yes	/	/	/
E7	SC	Yes	Yes	/	/	/
E8	SC	Yes	Yes	/	/	/
E9	SC	Yes	Yes	/	/	/
E10	SC	Yes	Yes	/	/	/
E11	SC	Yes	Yes	/	/	/
E12	SC	Yes	Yes	/	/	/
SD	SC	Yes	Yes	/	/	/
Number of participants		12				

10

NMR profiling
Quantitative parameters
1) Sugars & 5-HMF



- **Glucose + Fructose, sucrose, 5-HMF**
- **Comparison with directive 2001/110/EC**


- **Fructose / Glucose, Turanose**
- **Comparison with literature & databases**


AA0SG	AA	1H-NMR profiling of honey (sugars, HMF, sugar add)	Method : Internal, NMR	
(a)	Fructose	35.5 (± 2.8) g/100 g		
(a)	Glucose	31.0 (± 2.6) g/100 g		
(a)	Fructose / Glucose	1.15 (± 0.10)		
(a)	Glucose+Fructose	66.5 (± 4.3) g/100 g		>= 60 (blossom honey - miel de fleurs)
				>= 45 (Honeydew honey - miel de miellat) (Dir. 2001/110/CE)
				<= 5 (general ; Dir. 2001/110/CE)
(a)	Sucrose	<0.5 g/100 g		<= 40 (general ; Dir. 2001/110/CE)
(a)	5-HMF	22 (± 7) mg/kg		
(a)	Turanose	1.42 (± 0.19) g/100 g		

NMR profiling
Quantitative parameters
2) Others



- **citric acid**
- **fermentation markers (acetic, succinic & lactic acid, ethanol)**
- **Amino-acids**





- NMR profiling**
3) Detection of added sugars (both C3 & C4)
- **Analytical criteria for detecting sugar addition include:**
 - Concentration thresholds for typical honey compounds
 - Concentration thresholds for sugar syrup markers
 - Thresholds & Ratios between specific signals (ca 40 parameters)
 - Quantile plot for selected spectral windows
 - **Reporting:**
 - “The NMR profiling test”...
 - ...“did not reveal any peculiarities in this sample.”
 - or “shows an untypical spectrum, and points to an addition of exogenous sugar.”
 - **Interpretation harmonisation with Bruker consortium**
 - Information sharing
 - Case studies discussions
- 14


NMR profiling
4) Confirmation of the botanical origin



For monofloral honeys



Specific markers
 For botanical species



Complementary to pollen analysis
 (independent from potential pollen manipulation)

NMR profiling
3) Confirmation of the botanical origin
Manuka honey



Food Chemistry 217 (2017) 766–772

Contents lists available at ScienceDirect

Food Chemistry

journal homepage: www.elsevier.com/locate/foodchem




Analytical Methods

Combination of ¹H NMR and chemometrics to discriminate manuka honey from other floral honey types from Oceania

Marc Spiteri^a, Karyne M. Rogers^b, Eric Jamin^a, Freddy Thomas^{a,c}, Sophie Guyader^a, Michèle Lees^a, Douglas N. Rutledge^c


^aEurofins Analytics France, Rue Pierre Adolphe Bobierre, B.P. 42301, F-44323 NANTES Cedex 2, France
^bNational Isotope Center, CNS Science, 30 Grosvenor Road, Seaview, P.O. Box 31 312, Lower Hutt 5040, New Zealand
^cAgroParisTech, UMR1145 Ingénierie Procédés Aliments, 16, rue Claude Bernard, 75005 Paris, France



High value => fraud risk
Production << world consumption

- **Specific model to differentiate pure Manuka from non-100% manuka honey**
- **Specific test AA06C controls non peroxide activity**
- **High occurrence of non-compliances**

Stable isotope testing via EA-IRMS
Carbon 13: AOAC 998.12 official method



Difference of ^{13}C deviations




Detection and « quantification » of C4 plant (cane, maize) sugar
Limitation: does not detect C3 plants sugar (from wheat, rice, beet, etc.)

October 2008

17

Stable isotope testing via EA-IRMS
Carbon 13: LC-IRMS



Difference of ^{13}C deviations



Enhanced Detection of C4 plant (cane, maize) sugar
Limitation: does not detect C3 plants sugar (from wheat, rice, beet, etc.)

October 2008


18

Stable isotope testing via EA-IRMS
Geographical origin check
by multi-isotope analysis




eurofins

Isotopic fingerprint



Place of production



- Package PAA4W (¹³C, ¹⁵N, ²H), based on European database
- Test A4010 (⁸⁷Sr) for local origins based on geology data



Complementary to pollen analysis
(independent from potential pollen manipulation)


LC-HRMS
Advantages of High Resolution MS



eurofins


- Capable of separating mass fragments at the **fifth decimal place (exact mass)** where MS/MS instrumentation was limited to single-digit mass unit.
- high-resolution MS narrows down the possible chemical formulas of a unknown compound
- Example: Sucrose Mass
 $342.1 \rightarrow 274$ possible chemical formulas
 $342.11621 \rightarrow 1$ possible chemical formula
 $C_{12}H_{22}O_{11}$

Element	Nominal Mass	Accurate Mass
¹² C	12	12.00000000
¹ H	1	1.007825037
¹⁶ O	16	15.99491464
¹⁴ N	14	14.003074008





LC-HRMS

Specific Marker example : Difructose Anhydride




Tricyclic disaccharides formed during caramelization (14 known isomers)

HFCS and Inverted Syrup



Limit of Detection for Syrup Additions



**Case study of multi-methods approach:
Corsica honeydew honey**



- **Microscopic analysis**
(pollen, non soluble elements)
+Sensorial analysis



In agreement with the product description

- **C13-IRMS:**



No C4 sugar addition

- **NMR:**
very high citric acid content
(1327mg/kg)





Adulteration with sugar suspected

- **LC-HRMS:**
sugar syrup addition detected




Adulteration with sugar confirmed

**Non compliance statistics (all tests)
(From 2013 to 2019)**



All honey tests




Disclaimer :


Non compliance rates presented here are based on samples sent by customers and analyzed in Eurofins Analytics France. They may not reflect accurately the market situation in each country / segment

But they can be used as qualitative indications for major risks for each ingredient supply on a worldwide basis.

Non compliance statistics (all tests) (From 2013 to 2019)



NMR profiling only



Disclaimer :

Non compliance rates presented here are based on samples sent by customers and analyzed in Eurofins Analytics France. They may not reflect accurately the market situation in each country / segment

But they can be used as qualitative indications for major risks for each ingredient supply on a worldwide basis.

14/11/2019

Authenticity testing

25

Thank you for your attention !



Dr. Eric Jamin

Authenticity Competence Centre Manager

Eurofins Analytics France (Nantes)

14/11/2019

www.eurofins.com

www.eurofins.com