



## Food and Environmental Proficiency Testing Unit

The final date for return of results is indicated on the accompanying request/report form

### Contact details:

The Organisers - FEPTU  
Public Health England  
61 Colindale Avenue,  
London, NW9 5EQ, UK.  
Fax: +44 (0) 20 8200 8264  
Tel: +44 (0) 20 8327 7119  
e-mail: foodeqa@phe.gov.uk



0006

Results should be returned on-line using the web-based reporting system. Go to: [www.phe-qa.org.uk](http://www.phe-qa.org.uk). You will need your log-in details for this process. Please contact us if you experience problems with the on-line system.

[www.gov.uk/government/collections/external-quality-assessment-qa-and-proficiency-testing-pt-for-food-water-and-environmental-microbiology#legionella-isolation-scheme](http://www.gov.uk/government/collections/external-quality-assessment-qa-and-proficiency-testing-pt-for-food-water-and-environmental-microbiology#legionella-isolation-scheme)

## Water EQA Legionella Isolation – Instruction Sheet

Each distribution consists of **two** LENTICULE® discs in screw-cap plastic vials (with desiccant). The LENTICULE discs require reconstitution by a process of re-hydration and dispersion prior to examination for legionellae.

Refer also to the Safety Data Sheet: [www.gov.uk/government/publications/safety-data-sheet-for-lenticules](http://www.gov.uk/government/publications/safety-data-sheet-for-lenticules)

For a video on sample processing click on link below and scroll to bottom of page:  
[www.gov.uk/government/publications/legionella-scheme-sample-instruction-sheets](http://www.gov.uk/government/publications/legionella-scheme-sample-instruction-sheets)

### Storage:

- Store the samples at **-20 ± 5°C** on receipt.
- Allow the LENTICULE discs to reach ambient temperature (5 - 10 minutes) before reconstituting in diluent.

### Reconstitution:

- Open the sample container and transfer the LENTICULE disc into **1 Litre** of diluent (use either Page's saline or 1:40 Ringers' solution) by inverting the container over the diluent.
- Leave at ambient temperature for a minimum of 10 minutes to rehydrate. Ensure that the LENTICULE disc has dissolved completely before proceeding.
- Disperse the inoculum by inverting approximately 30 times.
- Examine samples between 30 - 45 minutes of reconstitution.

### Examination:

- Each reconstituted LENTICULE disc is equivalent to a 1 Litre water sample.**
- The request/report form indicates the sample dilutions most likely to give a definitive result (rather than a censored value).
- Once the dilutions have been prepared, the samples must be examined in accordance with routine procedures.
- All samples require examination for the presence of *Legionella* spp.
- The identification and level (colony forming units per Litre (cfu L<sup>-1</sup>)) of any legionellae detected must be determined.

### Return of Results:

- Return your results on-line using the web-based reporting system.
- The distribution will be closed on the deadline date for return of results.
- Results cannot be accepted after the distribution has closed.