


UK Standards for Microbiology Investigations

Identification of *Listeria* species, and other non-sporing gram positive rods (except *Corynebacterium*)



"NICE has renewed accreditation of the process used by Public Health England (PHE) to produce UK Standards for Microbiology Investigations. The renewed accreditation is valid until 30 June 2021 and applies to guidance produced using the processes described in UK standards for microbiology investigations (UKSMIs) Development process, S9365', 2016. The original accreditation term began in July 2011."

This publication was created by Public Health England (PHE) in partnership with the NHS.

Issued by the Standards Unit, National Infection Service, PHE

Bacteriology – Identification | ID 3 | Issue no: dm+ | Issue date: dd.mm.yy <tab+enter> |
Page: 1 of 21

Contents

Amendment table	3
1. General information.....	4
2. Scientific information	4
3. Scope of document	4
4. Introduction.....	4
5. Technical information/limitations.....	9
6. Safety considerations	9
7. Target organisms.....	10
8. Identification	10
9. Identification of <i>Listeria</i> species and other non-sporing Gram positive rods (except <i>Corynebacterium</i>).....	15
10. Referral to reference laboratories	16
11. Notification to PHE, or equivalent in the devolved administrations	17
References	18



"NICE has renewed accreditation of the process used by Public Health England (PHE) to produce UK Standards for Microbiology Investigations. The renewed accreditation is valid until 30 June 2021 and applies to guidance produced using the processes described in UK standards for microbiology investigations (UKSMIs) Development process, S9365', 2016. The original accreditation term began in July 2011."

Amendment table

Each UK SMI method has an individual record of amendments. The current amendments are listed on this page. The amendment history is available from standards@phe.gov.uk.

New or revised documents should be controlled within the laboratory in accordance with the local quality management system.

Amendment number/date	New amendment number/dd.mm.yy <tab+enter>
Issue number discarded	
Insert issue number	
Anticipated next review date*	
Section(s) involved	Amendment
Whole document.	Taxonomy of <i>Listeria</i> species updated. Minor amendments to characteristics of organisms. Information on whole genome sequencing added.

*Reviews can be extended up to five years subject to resources available.

1. General information

[View](#) general information related to UK SMIs.

2. Scientific information

[View](#) scientific information related to UK SMIs.

3. Scope of document

This UK Standard for Microbiology Investigation (UK SMI) describes the identification of *Listeria* species and other non-sporing Gram positive rods (except *Corynebacterium* species) isolated from clinical specimens to genus or species level.

For the identification of *Corynebacterium* species, refer to [ID 2 - Identification of *Corynebacterium* species](#). This UK SMI should be used in conjunction with other UK SMIs.

4. Introduction

A systematic approach is used to differentiate clinically encountered, morphologically similar, aerobic and facultatively anaerobic, non-sporing Gram positive rods. The true branching organisms such as *Actinomyces*, *Nocardia* and *Streptomyces* species and those which produce spores are not described in this UK SMI. Rapidly growing *Mycobacterium* species may also be isolated on the media described in this document and acid-fast bacilli should be referred to the Reference Laboratory.

4.1 Taxonomy/characteristics

Listeria species¹⁻⁴

There are currently nineteen validly named species in the genus *Listeria*:

L. monocytogenes, *L. ivanovii*, *L. seeligeri*, *L. innocua*, *L. welshimeri*, *L. grayi*, *L. fleischmannii*, *L. marthii*, *L. rocourtiae*, *L. weihenstephanensis*, *L. aquatica*, *L. booriae*, *L. cornellensis*, *L. floridensis*, *L. grandensis*, *L. newyorkensis*, *L. costaricensis*, *L. murrayi* and *L. riparia* of these nineteen species, the first six can potentially cause infections in humans, albeit rarely in some cases.

The type species is *Listeria monocytogenes*.

Listeria species are short Gram positive rods, 0.4-0.5 x 0.5-2.0 µm, with rounded ends, occurring singly or in short chains. Filaments of 6 to 20 µm in length may occur singly or in short chains. Members of the genus are facultative anaerobes, non-sporing, non-acid fast and do not possess a capsule. The optimum growth temperature (but not for motility) is 30-37°C.

Colonies on blood agar are non-pigmented and may resemble those of β-haemolytic streptococci but can be distinguished by a positive catalase test. Colonies will vary in colour and form depending on which commercial chromogenic agar is being used.

Listeria species are oxidase negative and ferment carbohydrates, they are widely distributed in the environment; some species are pathogenic for humans and animals.

The medically important species are:

***L. monocytogenes*³**

Microscopically, they appear as small rods, which are sometimes arranged in short chains. In direct smears they may be coccoid, and may be mistaken for streptococci. Longer rods may resemble *Corynebacterium*. Haemolytic activity on blood agar has been used as a marker to distinguish *Listeria monocytogenes* from other *Listeria* species, but it is not an absolutely definitive criterion. Further biochemical characterization may be necessary to distinguish between the different *Listeria* species. They are oxidase negative.

L. monocytogenes is the agent of listeriosis, a serious infection caused by eating food contaminated with the bacterium. Listeriosis has been recognized as an important public health problem and the disease affects primarily pregnant women, neonates, elderly people, and those with weakened immune systems.

***L. ivanovii*⁵**

This species has been divided into 2 subspecies. These are; *Listeria ivanovii* subsp. *ivanovii* and *Listeria ivanovii* subsp. *Londoniensis*⁶. They are facultatively anaerobic and has also been isolated from healthy animals and human carriers and from the environment.

***L. seeligeri*⁷**

Colonies on tryptose agar are similar to that of *L. welshimeri*. Growth occurs at 4°C within 5 days. They are facultatively anaerobic.

They have been isolated from plants, soil, and animal faeces (sheep) in Europe.

***L. innocua*^{8,9}**

They are mesophilic, operating at an optimal temperature range of 30-37°C. *Listeria innocua* have a very complex metabolism. They are capable of metabolizing methane, sulphur and nitrogen, among many other organic and inorganic compounds. These organisms also carry out numerous biosynthetic pathways, including peptidoglycan synthesis. *L. innocua*, like other members of their genus, are facultative anaerobes, which means that they can metabolism glucose (and other simple sugars) in under both aerobic and anaerobic conditions. Under the aerobic metabolism of glucose, they form lactic acid and acetic acid. However, under anaerobic conditions, the metabolism of glucose yields only lactic acid.

This species is widespread in the environment and in food and has also been associated with one reported case of fatal bacteraemia.

***L. welshimeri*⁷**

Cells are small (0.4 to 0.5 by 0.5 to 2.0µm) rods which are motile by means of peritrichous flagella. Colonies on tryptose agar are small (1 to 2mm in diameter after 1 or 2 days of incubation at 37°C), regular, and smooth with a blue-green colour when they are examined with obliquely transmitted light. Sheep erythrocytes are not haemolysed. Growth occurs at 4°C within 5 days.

Metabolism is facultatively anaerobic. Acid, but no gas, is produced from D-glucose, D-xylose, and α-methyl-D-mannoside. Acid may or may not be produced from L- rhamnose. Acid is not produced from D-mannitol. They are, aesculin hydrolysis,

Voges-Proskauer and methyl red tests, and negative for oxidase, urea, gelatin hydrolysis, indole and H₂S production as well as reduction of nitrates.

They have been isolated from decaying plants and soil.

***L. grayi*¹⁰**

According to Rocourt *et al.* (1992), *Listeria grayi* is an earlier heterotypic synonym of *Listeria murrayi* and so both were assigned to a single species, *Listeria grayi*.

L. grayi peritrichous rods which are motile. Colonies on tryptose agar are small and growth occurs at 4°C within 5 days.

Metabolism is facultatively anaerobic. They are, aesculin hydrolysis, Voges-Proskauer and methyl red tests and negative for the oxidase, urea and gelatin hydrolysis, H₂S and indole production. Reduction of nitrates to nitrites is variable. Acid, without gas, is produced from glucose, mannitol, and other sugars. Sheep erythrocytes are not haemolysed.

Other non-sporing Gram positive rods¹¹⁻¹³

The organisms classified as non-sporing Gram positive rods are very diverse not only morphologically, but also metabolically and structurally.

***Arcanobacterium* species^{14,15}**

There were 13 species of which 5 have been re-assigned to the genus *Trueperella* and of the remaining 8 species, only one is known to infect humans, *Arcanobacterium haemolyticum*.

Both rod-shaped and coccoid cells are Gram positive, non-acid-fast and non-motile. Endospores are not formed. They are facultatively anaerobic. Growth is considerably enhanced in an atmosphere of CO₂. Growth is sparse on ordinary media but enhanced on blood or serum containing media. Optimum temperature for growth is 37°C. They are unable to withstand heating at 60°C for 15 min.

Arcanobacterium species may give variable catalase reactions and are positive for CAMP-test. See species below.

***Arcanobacterium haemolyticum* (formerly *Corynebacterium haemolyticum*)¹⁶**

Colonies on blood agar after 48hr produce zones of β-haemolysis and are similar in appearance to *Trueperella pyogenes*.

They have been isolated from the throat of infected individuals.

***Aureobacterium* species^{17,18}**

Aureobacterium species are obligate aerobes, which produce acid from carbohydrates by oxidation rather than by fermentation. Strains may be vancomycin resistant and can be distinguished from *C. aquaticum* by casein and gelatin hydrolysis.

***Bifidobacterium* species¹**

Bifidobacterium species vary in shape and may be curved, clubbed or branched rods or occasionally coccoid, Gram positive forms, 0.5 - 1.3 x 1.5-8µm. Cells often stain irregularly. Growth is anaerobic but some species can grow in air enriched with 10% CO₂. *Bifidobacterium* species are non-sporing, non-acid fast and ferment carbohydrates.

***Brevibacterium* species^{16,19}**

Brevibacterium species are non-motile, salt tolerant (>6.5% NaCl), aerobic, urease negative and catalase positive.

***Cellulomonas* species^{16,20}**

There are over 20 *Cellulomonas* species, of which only *C. hominis* and *C. denverensis* have been isolated from humans. They may be non-motile or motile due to single or sparse lateral flagella. One of their main distinguishing features is their ability to degrade cellulose, using enzymes such as endoglucanase and exoglucanase. They are both oxidative and fermentative in their metabolism. *Cellulomonas* species differ from *Oerskovia* species in that they lack hyphal growth.

***Cutibacterium* (formerly *Propionibacterium* species)²¹**

The genus *Cutibacterium* has only 3 species which have been re-assigned from the genus *Propionibacterium* and they are *C. acnes*, *C. avidum* and *C. granulosum*. They are Gram positive pleomorphic non-motile rods and generally grow better anaerobically, producing small colonies after 5 - 6 days.

***Dermabacter hominis*¹⁶**

Dermabacter hominis, currently the only member of the genus, is non-haemolytic, non-motile and catalase positive. *Dermabacter* species are fermentative and produce acid from glucose, lactose, sucrose and maltose. They hydrolyse aesculin and produce alkaline phosphatase, pyrrolidonyl arylamidase, leucine aminopeptidase and DNase. They do not reduce nitrate or produce pyrazinamidase.

***Erysipelothrix rhusiopathiae*²²**

E. rhusiopathiae produces a narrow zone of α -haemolysis on blood agar. It is facultatively anaerobic, non-motile and catalase negative. All colonies are clear, circular and very small increasing in size and tending towards a pale blue opacity with further incubation or age. *Erysipelothrix* species can be distinguished from *Lactobacillus* species by its ability to produce H₂S in a triple sugar iron agar slant.

***Gardnerella vaginalis*¹**

Gardnerella vaginalis is a pleomorphic, Gram variable rod or coccobacilli. It is facultatively anaerobic and non-motile. *G. vaginalis* is non-sporing, non-encapsulated and both oxidase negative and catalase negative. Acid is produced from glucose and other carbohydrates but not gas. It hydrolyses hippurate and does not reduce nitrate.

***Lactobacillus* species¹**

Lactobacillus species are long Gram positive rods. Colonies can vary from small to medium, grey and often α -haemolytic on blood agar after 48hr. They are facultatively anaerobic, rarely motile and catalase negative. They produce lactic acid as a major end product.

***Microbacterium* species^{16,20}**

Microbacterium species are small, slender, irregularly shaped Gram positive rods. The optimum growth temperature is 30°C. The species are primarily oxidative and aerobic in their metabolism, but some species may be fermentative. They may be non-motile or motile by means of 1 to 3 flagella. Most species are catalase positive but catalase negative strains may be observed.

***Mycobacterium* species**

Mycobacterium species other than *Mycobacterium tuberculosis* (MOTT) may be isolated on primary culture within 48hr for identification and/or susceptibility ([B 40 - Investigation of specimens for *Mycobacterium* species](#)). Refer to the National Mycobacterium Reference Laboratory.

***Oerskovia* species¹⁶**

Oerskovia species are Gram positive branching rods. They form a mycelium, an extensively branching substrate hypha that breaks up to form rod-shaped motile or non-motile coccoid-rod elements. No aerial hyphae are formed. They are facultatively anaerobic, fermentative and catalase positive. Motility is variable.

***Propionibacterium* species¹⁶**

Propionibacterium species are Gram positive pleomorphic rods (short “Y” forms). Strains generally grow better anaerobically, particularly on primary isolation, producing small colonies after 48hr. *Propionibacterium* species are facultatively anaerobic and are non-motile. They are catalase positive except *Propionibacterium propionicum* (now known as *Pseudopropionibacterium propionicum*), which is catalase negative.

Note: It should be noted that within this genus, some of the species have been reclassified to other genera.

***Rhodococcus* species²³**

Rhodococcus species usually stain Gram positive. Cells form as cocci or short rods which grow in length, and may form an extensively branched vegetative mycelium which may fragment. They are usually partially acid-fast due to the mycolic acid in their cell walls. Incubation at 30°C also increases recovery.

Although other biochemical tests help to distinguish *Rhodococcus* from other organisms, differentiation from other aerobic actinomycetes can be difficult. *Rhodococcus* species typically react positively in catalase, nitrate reduction, and urea hydrolysis tests and negatively with oxidase, gelatin hydrolysis, and carbohydrate reduction. They are non-motile and non-sporulating. Their inability to ferment carbohydrate is important in distinguishing them from corynebacteria.

***Trueperella* species^{14,24,25}**

There are currently five validly named species in the genus *Trueperella* and two of these cause infections in humans.

Cells are, non-motile, non-spore-forming coccobacilli and. Cells vary in shape and size (0.2–0.9 x 60.3–2.5µm) in different media. Cells from 24hour old broth cultures are Gram positive, but may be Gram variable in older cultures. Pinpoint, β-haemolytic colonies occur on sheep's blood agar after 24hr of incubation. They are aerobic and facultatively anaerobic. Members are strictly fermentative. Lactic acid is the primary metabolic product in glucose yeast extract broth but acetate and succinate are minor products.

The type species is *Trueperella pyogenes*.

***Trueperella pyogenes*²⁶** (formerly *Arcanobacterium pyogenes*) is a rod which may show branching.

T. pyogenes is facultatively anaerobic, non-motile, and catalase negative but one strain has been reported as positive. Metabolism is strictly fermentative. Differentiation between *T. pyogenes* and *A. haemolyticum* may prove difficult but they may be distinguished by fermentation of α -mannose, pyrazinamidase and gelatin tests.

This organism is frequently isolated from a wide variety of pyogenic disease conditions in many species of domestic animals and in humans.

***Trueperella bernardiae*²⁷** (formerly *Arcanobacterium bernardiae*) is a rod with coccobacilli predominating. Primary branching is not observed. *T. bernardiae* is facultatively anaerobic, non-motile, catalase negative and unable to produce acid from glucose. after 48hr incubation. Colony diameters range from 0.2mm to 0.5mm after 48hr incubation.

This has been isolated from human blood as well as abscess from ear or chest.

***Turicella otitidis*^{16,20}**

The genus comprises a single species, *Turicella otitidis*. Microscopically it resembles a coryneform but has longer cells. It may be distinguished by colonial morphology from *Corynebacterium afermentans* and *Corynebacterium auris*. compared with the flat, grey-white and non-haemolytic colonies of *C. afermentans* and the convex, dry, adherent, yellowish colonies of *C. auris*. *T. otitidis* is non-fermentative and occurs either alone or with Gram negative rods. Isolates exhibit a strong CAMP reaction and are DNase positive and catalase positive.

T. otitidis may be misidentified, often as *Corynebacterium* species, by some commercial identification systems.

5. Technical information/limitations

Differentiation of *Listeria* species from Group B streptococci

Colonies of *Listeria* species resemble those of Group B streptococci, and the catalase test is a rapid, easily performed test which will help differentiate *Listeria* species from Group B streptococci. *Listeria* species are catalase positive whereas Group B streptococci are catalase negative.

Motility test (see below)

Motility is one of many parameters used in the characterisation of *Listeria* species. It should be used in conjunction with other tests. This test should not be used for primary isolation of *Listeria* species or purposes other than the investigation of motility.

Listeria species are motile by peritrichous flagella when at 20-28°C and non-motile at 35°C and above. Therefore, an appropriate temperature must be chosen for incubation to avoid false negative results. There have been occasional non-motile strains²⁸.

6. Safety considerations^{11-13,29-43}

Hazard Group 2 organisms

Pregnant staff should be prohibited from working with known or suspected cultures of *Listeria* species.

Laboratory procedures that give rise to infectious aerosols must be conducted in a microbiological safety cabinet¹³.

Refer to current guidance on the safe handling of all organisms documented in this UK SMI.

The above guidance should be supplemented with local COSHH and risk assessments.

Compliance with postal and transport regulations is essential.

7. Target organisms

***Listeria* species and other morphologically similar Gram positive rods reported to have caused human infection^{16,44,45}**

Listeria monocytogenes, *Arcanobacterium* species, *Erysipelothrix rhusiopathiae*

Other *Listeria* species reported to have caused occasional or single human infections^{9,46-51}

Listeria ivanovii, *Listeria seeligeri*, *Listeria innocua*, *Listeria grayi*

Other species morphologically similar to *Listeria* species known to have caused human infection^{11-13,16}

Aureobacterium species, *Bifidobacterium* species, *Brevibacterium* species, *Cellulomonas* species, *Dermabacter hominis*

Other Gram positive rods have been implicated in human infections^{16,20,52}.

8. Identification

8.1 Microscopic appearance

Gram stain ([TP 39 - Staining procedures](#))

Gram positive rods. Microscopic appearance varies with the species.

***Listeria* species**

Gram positive rods approximately 0.5 x 0.5 – 3µm with rounded ends, occurring singly or sometimes in pairs and may resemble 'coryneforms' or diplococci. They are non-sporing, non-branching and non-capsulated.

***Arcanobacterium* species**

Gram positive rod-shaped and coccoid cells; they are slender, irregular and predominately rod-shaped or arranged at an angle to give V-formations during the first 18hr of growth, becoming granular and segmented, resembling small, irregular cocci over time.

Erysipelothrix rhusiopathiae

Cells are slender Gram positive non-sporulating rods occurring in short chains, in pairs, in a "V" configuration or even grouped randomly. This organism can appear Gram negative because of their tendency to decolourise rapidly.

***Trueperella* species**

Cells are Gram positive, coccobacilli and short rods that occur singly, in pairs (V, T and palisade formations) or in clusters.

***Aureobacterium* species**

Cells are Gram positive irregular short rods.

***Bifidobacterium* species**

These species are Gram positive, vary in shape and may be curved, clubbed or branched rods or occasionally coccoid forms, 0.5-1.3 x 1.5-8µm. Cells often stain irregularly.

***Brevibacterium* species**

Gram positive rods and they show a marked rod-coccus cycle. On fresh subculture, cells appear as bacilli but become coccal in older cultures.

***Cellulomonas* species**

Cells are Gram positive slender irregular rods that produce yellow or orange pigmented colonies

Dermabacter hominis

They are very short Gram positive rods that may be misinterpreted as cocci.

For all the other non-sporing Gram positive rods, see the section for “Characteristics”.

8.2 Primary isolation media

Blood agar incubated in 5–10% CO₂ at 35°C–37°C for 16–48hr.

Listeria selective agar incubated in O₂ at 35°C–37°C for 40–48hr.

Note: *Listeria* species are also capable of growth at 2°C–43°C.

8.3 Colonial appearance

This table is a summary of the organisms and their appearances on blood agar plate.

Organism	Characteristics of growth on blood agar after incubation at 35-37°C for 16-48hr
<i>L. monocytogenes</i>	Colonies are 0.5-1.5mm in diameter, smooth, translucent with a characteristic ground glass appearance able to be emulsified and with a zone of hazy β-haemolysis extending 1-2mm from the edge of the colony.
<i>Arcanobacterium</i> species	Cells are slender, irregular and predominately rod-shaped or arranged at an angle to give V-formations during the first 18hr of growth on blood agar, becoming granular and segmented, resembling small, irregular cocci over time. After 48hr incubation, colonies produce zones of β-haemolysis.
<i>Erysipelothrix rhusiopathiae</i>	After 48hr incubation, two distinct colony types appear: a small smooth (S) form, 0.5-1mm in diameter, transparent, convex and circular with entire edges. The large rough (R) form is flatter, more opaque, with a matt surface and an irregular edge. Most strains exhibit a narrow zone of α-haemolysis but the R –form does not cause haemolysis.

<i>L. ivanovii</i>	Colonies are similar to <i>L. monocytogenes</i> but develop larger zones of complete haemolysis with outer zones of partial haemolysis. Cells are small, motile rods. Colonies on tryptose agar are very small (0.5 to 1mm in diameter after 1 or 2 days of incubation at 37°C), regular, and smooth and appear bluish green when they are viewed by obliquely transmitted light. Colonies on sheep or horse blood (5%) agar are strongly β-haemolytic. Growth occurs at 4°C within 5 days.
<i>L. seeligeri</i>	Cells are small (0.4 to 0.8 x 0.5 to 2.5µm) rods Colonies are similar to <i>L. monocytogenes</i> but zones of β- haemolysis are produced.
<i>L. innocua</i>	Cells are small rods occurring singly or in short chains. Cream colonies, no haemolysis.
<i>L. grayi</i>	Colonies are small, regular, smooth and are 1 to 2mm in diameter after 1- 2 days of incubation at 37°C.
<i>Trueperella</i> species	Pinpoint, β-haemolytic colonies occur on blood agar after 24hr of incubation and after 48–72hr of incubation, colonies are 0.5–1.5mm in diameter, convex, circular and translucent with entire edges. <i>T. pyogenes</i> produces sharp zones of β-haemolysis after 48hr incubation. Haemolysis of <i>T. bernardiae</i> is variable.
<i>Aureobacterium</i> species	Non-haemolytic, yellow pigmented colonies.
<i>Bifidobacterium</i> species	Colonies are low, greyish-brown, and ovoid with a brown opaque centre and translucent crenated edges.
<i>Brevibacterium</i> species	Colonies are opaque, grey-white, 2mm or more in diameter after 24hr, convex and have a smooth shiny surface. They are non-haemolytic and may turn yellow to green after 48hr.
<i>Cellulomonas</i> species	Non-haemolytic, yellow- or orange-pigmented colonies.
<i>Dermabacter hominis</i>	Non-haemolytic, small grey/white convex colonies with entire edges.
<i>Lactobacillus</i> species	Colonies are small and often α-haemolytic on blood agar after 48hr.
<i>Microbacterium</i> species	Colonies are circular, convex with entire margins, moist, shiny and may produce a yellow or orange pigment.
<i>Oerskovia</i> species	Most strains produce a yellow pigment.
<i>Propionibacterium</i> species	They produce small colonies after 48hr incubation.
<i>Gardnerella vaginalis</i>	Growth is enhanced by the addition of 5-10% CO ₂ Colonies are small, circular, convex and grey. It also produces diffuse β-haemolysis on rabbit blood agar but not on sheep blood agar. Haemolysis on horse blood agar is variable.
<i>Rhodococcus</i> species	Colonies may be rough, smooth or mucoid and are colourless, cream, beige, yellow, orange or red.
<i>Turicella otitidis</i>	Colonies are convex, whitish, creamy and non-haemolytic in appearance.
<i>Trueperella bernardiae</i>	Colonies exhibit variable haemolysis, appear circular, smooth and slightly convex with a glassy appearance.
<i>Trueperella pyogenes</i>	Colonies produce sharp zones of β-haemolysis after 48hr incubation. They appear convex, white, smooth, translucent and soft with entire edges.

Other *Listeria* species show similar colonial appearance and are either haemolytic or non-haemolytic: these species are very rarely isolated from normally sterile sites and should be submitted to the Reference Laboratory for identification.

For all other non-sporing Gram positive rods^{16,20}

Appearance varies with species on blood agar, after aerobic incubation at 35–37°C for 16–48hr. See table above.

8.4 Test procedures

Catalase test ([TP 8 - Catalase test](#))

Listeria species are catalase positive.

Arcanobacterium species may give variable catalase reactions.

Erysipelothrix rhusiopathiae is catalase negative.

For the other non-sporing rods, see the flowchart and table 1 on section 8.

Motility test ([TP 21 - Motility test](#))

This is performed at 20°C - 28°C for *Listeria* species and above 30°C for all other organisms.

All *Listeria* species exhibit tumbling motility at 20°C - 28°C but not at above 30°C. Other organisms may be motile but do not exhibit tumbling motility.

For the other non-sporing rods, see the flowchart and table 1 on section 8.

Commercial identification systems⁵³

Laboratories should follow manufacturer's instructions. Rapid tests and kits and should be validated and be shown to be fit for purpose prior to use.

Matrix Assisted Laser Desorption Ionisation Time-of-Flight (MALDI-TOF)⁵⁴

This has been shown to be a rapid and powerful tool because of its reproducibility, speed and sensitivity of analysis.

MALDI-TOF has been developed and validated to determine species and lineages of *Listeria* species using isolates from a variety of sources.

It has been known to be used to identify *T. bernardiae* and thus will help in its future identification and in elucidating the role that this rarely isolated species plays in infection of humans⁵⁵.

Refer to UK SMI [TP 40 - MALDI TOF MS test procedure](#) for further information.

Nucleic Acid Amplification Tests (NAATs)^{56,57}

PCR is usually considered to be a good method as it is simple, sensitive and specific. The basis for PCR diagnostic applications in microbiology is the detection of infectious agents and the discrimination of non-pathogenic from pathogenic strains by virtue of specific genes.

This is now established as a rapid, reliable and reproducible technique for identification of *Listeria* species and, more importantly, for the differentiation of

L. monocytogenes from other species using primers targeting genes encoding virulence factors or RNA sub-unit genes.

The Multiplex PCR-based serotyping assay, such as the use of group-specific PCR primers, has equally provided additional tools for the identification and grouping of *L. monocytogenes*.

8.5 Further identification

Following the colonial morphology, catalase test, motility test and biochemical identification results, if further identification is required, send isolate to the Reference Laboratory.

Whole genome sequencing (WGS)⁵⁸

All isolates referred to the reference laboratory are subjected to whole genome sequencing (WGS) also known as full genome sequencing, complete genome sequencing or entire genome sequencing is a laboratory method that determine the complete DNA sequence of an organism's genome at a single time.

WGS has been used successfully to explore the genome of organisms as well as to identify candidate genes responsible for pathogenesis, and to develop better methods of strain identification and to advance the understanding of evolution of organisms. WGS analysis can provide phylogenetic data, which helps to perform real-time or retrospective surveillance of *L. monocytogenes*.

Single nucleotide polymorphism (SNP)⁵⁹


The single nucleotide polymorphism (SNP) based approach is most commonly used principle for the retrieval of subtyping information from WGS data. SNPs are flexible, do not require a redefined scheme and provides an exceptionally high subtyping resolution.

SNP is used for the highly discriminative typing of *L. monocytogenes*.

8.6 Storage and referral

If required, save the pure isolate on a blood or nutrient agar slope for referral to the Reference Laboratory.

9. Identification of *Listeria* species and other non-sporing Gram positive rods (except *Corynebacterium*)



The flowchart is for guidance only.

Table 1: Phenotypic characteristics of *Listeria* and other non-sporing gram positive rods (except *Corynebacterium*)

Organisms	Gram stain	Catalase	Motility
<i>Listeria</i> Species	Positive	Positive	Positive
<i>Arcanobacterium haemolyticum</i>	Positive	Negative	Negative
<i>Aureobacterium</i> species	Positive	Positive	Variable
<i>Bifidobacterium</i> species	Positive	Negative	Negative
<i>Brevibacterium</i> species	Positive	Positive	Negative
<i>Cellulomonas</i> species	Positive	Positive	Variable
<i>Dermabacter hominis</i>	Positive	Positive	Negative
<i>Erysipelothrix rhusiopathiae</i>	Positive	Negative	Negative
<i>Gardnerella vaginalis</i>	Variable	Negative	Negative
<i>Lactobacillus</i> species	Positive	Negative	Rarely motile
<i>Microbacterium</i> species	Positive	Variable	Variable
<i>Oerskovia</i> species	Positive	Positive	Variable
<i>Propionibacterium</i> species	Positive	Positive	Negative
<i>Rhodococcus</i> species	Positive	Positive	Negative
<i>Trueperella pyogenes</i>	Positive	Negative	Negative
<i>Trueperella bernardiae</i>	Positive	Negative	Negative
<i>Turicella otitidis</i>	Positive	Positive	Negative

10. Referral to reference laboratories

Report susceptibilities as clinically indicated. Prudent use of antimicrobials according to local and national protocols is recommended.

For information on the tests offered, turnaround times, transport procedure and the other requirements of the reference laboratory [click here for user manuals and request forms](#).

Organisms with unusual or unexpected resistance, or associated with a laboratory or clinical problem, or anomaly that requires elucidation should be sent to the appropriate reference laboratory.

Contact appropriate devolved national reference laboratory for information on the tests available, turnaround times, transport procedure and any other requirements for sample submission:

England and Wales

<https://www.gov.uk/specialist-and-reference-microbiology-laboratory-tests-and-services>

Scotland

<http://www.hps.scot.nhs.uk/reflab/index.aspx>

Northern Ireland

<http://www.publichealth.hscni.net/directorate-public-health/health-protection>

11. Notification to PHE^{60,61}, or equivalent in the devolved administrations⁶²⁻⁶⁵

The Health Protection (Notification) regulations 2010 require diagnostic laboratories to notify Public Health England (PHE) when they identify the causative agents that are listed in Schedule 2 of the Regulations. Notifications must be provided in writing, on paper or electronically, within seven days. Urgent cases should be notified orally and as soon as possible, recommended within 24 hours. These should be followed up by written notification within seven days. All other reports should be made as soon as possible after the causative agent has been identified, preferably at the same time that the laboratory report is being issued to the requester.

For the purposes of the Notification Regulations, the recipient of laboratory notifications is the local PHE Health Protection Team. If a case has already been notified by a registered medical practitioner, the diagnostic laboratory is still required to notify the case if they identify any evidence of an infection caused by a notifiable causative agent.

Notification under the Health Protection (Notification) Regulations 2010 does not replace voluntary reporting to PHE. The vast majority of NHS laboratories voluntarily report a wide range of laboratory diagnoses of causative agents to PHE and many PHE Health protection Teams have agreements with local laboratories for urgent reporting of some infections. This should continue.

Note: The Health Protection Legislation Guidance (2010) includes reporting of Human Immunodeficiency Virus (HIV) & Sexually Transmitted Infections (STIs), Healthcare Associated Infections (HCIs) and Creutzfeldt–Jakob disease (CJD) under ‘Notification Duties of Registered Medical Practitioners’: it is not noted under ‘Notification Duties of Diagnostic Laboratories’.

<https://www.gov.uk/government/organisations/public-health-england/about/our-governance#health-protection-regulations-2010>

Other arrangements exist in [Scotland](#)^{62,63}, [Wales](#)⁶⁴ and [Northern Ireland](#)⁶⁵.

References

1. Holt JG. Regular, nonsporing gram-positive rods. Genus *Listeria*. In: Holt JG, Krieg NR, Sneath PHA, Staley JT, Williams ST, editors. *Bergey's Manual of Determinative Bacteriology*. 9th ed. Baltimore: Williams and Wilkins; 1994. p. 565-70.
2. Euzéby JP. Genus *Listeria*. 2013.
3. Patrick R. Murray EJB, James Jorgensen, Michael Pfaller, Marie Louise Landry. *Manual of Clinical Microbiology* 2007;9.
4. Orsi RH, Wiedmann M. Characteristics and distribution of *Listeria* spp., including *Listeria* species newly described since 2009. *Applied microbiology and biotechnology* 2016;100:5273-87.
5. Seeliger HPR, Rocourt J, Schrettenbrunner A, Grimont PA, D J. *Listeria ivanovii* sp. nov. *International Journal of Systematic Bacteriology* 1984;34:336-7.
6. Boerlin P, Rocourt J, Grimont F, Grimont PA, Jacquet C, Piffaretti JC. *Listeria ivanovii* subsp. *londoniensis* subsp.nov. *International Journal of Systematic Bacteriology* 1992;42:69-73.
7. Rocourt J, Grimont PA. *Listeria welshimeri* sp. nov. and *Listeria seeligeri* sp.nov. *International Journal of Systematic Bacteriology* 1983;33:866-9.
8. Pine L, Malcolm GB, Brooks JB, Daneshvar MI. Physiological studies on the growth and utilization of sugars by *Listeria* species. *CanJMicrobiol* 1989;35:245-54.
9. Perrin M, Bemer M, Delamare C. Fatal case of *Listeria innocua* bacteremia. *JClinMicrobiol* 2003;41:5308-9.
10. Rocourt J, Boerlin P, Grimont F, Jacquet C, Piffaretti JC. Assignment of *Listeria grayi* and *Listeria murrayi* to a single species, *Listeria grayi*, with a revised description of *Listeria grayi*. *IntJSystBacteriol* 1992;42:171-4.
11. Advisory Committee on Dangerous Pathogens. The Approved List of Biological Agents. Health and Safety Executive 2013. 1-35. **A, VI**
12. Advisory Committee on Dangerous Pathogens. Infections at work: Controlling the risks. Her Majesty's Stationery Office 2003. **A, VI**
13. Advisory Committee on Dangerous Pathogens. Biological agents: Managing the risks in laboratories and healthcare premises. Health and Safety Executive 2005. **A, VI**
14. Yassin AF, Hupfer H, Siering C, Schumann P. Comparative chemotaxonomic and phylogenetic studies on the genus *Arcanobacterium* Collins et al. 1982 emend. Lehtinen et al. 2006: proposal for *Trueperella* gen. nov. and emended description of the genus *Arcanobacterium*. *IntJSystEvolMicrobiol* 2011;61:1265-74.
15. Euzéby JP. Genus *Arcanobacterium*. 2013.
16. Funke G, von Graevenitz A, Clarridge JE, III, Bernard KA. Clinical microbiology of coryneform bacteria. *Clin Microbiol Rev* 1997;10:125-59.

17. Nolte FS, Arnold KE, Sweat H, Winton EF, Funke G. Vancomycin-resistant *Aureobacterium* species cellulitis and bacteremia in a patient with acute myelogenous leukemia. *JClinMicrobiol* 1996;34:1992-4.
18. Funke G, von GA, Weiss N. Primary identification of *Aureobacterium* spp. isolated from clinical specimens as "*Corynebacterium aquaticum*". *JClinMicrobiol* 1994;32:2686-91.
19. Gruner E, Pfyffer GE, von GA. Characterization of *Brevibacterium* spp. from clinical specimens. *JClinMicrobiol* 1993;31:1408-12.
20. Janda WM. *Corynebacterium* species and the Coryneform bacteria part II: current status of the cdc coryneform groups. *Clin Microbiol Newslett* 1998;20:53-66.
21. Scholz CF, Kilian M. The natural history of cutaneous propionibacteria, and reclassification of selected species within the genus *Propionibacterium* to the proposed novel genera *Acidipropionibacterium* gen. nov., *Cutibacterium* gen. nov. and *Pseudopropionibacterium* gen. nov. *International journal of systematic and evolutionary microbiology* 2016;66:4422-32.
22. Brooke CJ, Riley TV. *Erysipelothrix rhusiopathiae*: bacteriology, epidemiology and clinical manifestations of an occupational pathogen. *JMed Microbiol* 1999;48:789-99.
23. Conlin M, Willingham K, Oliver J. *Rhodococcus* Bacteremia in an Immunocompetent Patient. *Laboratory Medicine* 2000;31:263-5.
24. Ramos CP, Foster G, Collins MD. Phylogenetic analysis of the genus *Actinomyces* based on 16S rRNA gene sequences: description of *Arcanobacterium phocae* sp. nov., *Arcanobacterium bernardiae* comb. nov., and *Arcanobacterium pyogenes* comb. nov. *IntJSystBacteriol* 1997;47:46-53.
25. Euzéby JP. Genus *Trueperella*. 2013.
26. Reddy CA, Cornell CP, Fraga AM. Transfer of *Corynebacterium pyogenes* (Glage) Ebersson to the Genus *Actinomyces* as *Actinomyces pyogenes* (Glage) comb. nov. *International Journal of Systematic Bacteriology* 1982;32:419-29.
27. Funke G, Ramos CP, Fernandez-Garayzabal JF, Weiss N, Collins MD. Description of human-derived Centers for Disease Control coryneform group 2 bacteria as *Actinomyces bernardiae* sp. nov. *IntJSystBacteriol* 1995;45:57-60.
28. Mitchell RG. *Listeria* and *Erysipelothrix*. In: Collee JG, Fraser AG, Marmion BP, Simmons A, editors. *Mackie & McCartney Practical Medical Microbiology*. 14th ed. Edinburgh: Churchill Livingstone; 1996. p. 309-15.
29. European Parliament. UK Standards for Microbiology Investigations (UK SMIs) use the term "CE marked leak proof container" to describe containers bearing the CE marking used for the collection and transport of clinical specimens. The requirements for specimen containers are given in the EU *in vitro* Diagnostic Medical Devices Directive (98/79/EC Annex 1 B 2.1) which states: "The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes". 1998. **A, VI**
30. Official Journal of the European Communities. Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices 1998. 1-37. **A, VI**

31. Health and Safety Executive. Safe use of pneumatic air tube transport systems for pathology specimens. 2009. **A, VI**
32. Department for Transport. Transport of Infectious Substances, 2011 Revision 5. 2011. **A, VI**
33. World Health Organization. Guidance on regulations for the Transport of Infectious Substances 2017-2018. 2017. **A, VI**
34. Home Office. Anti-terrorism, Crime and Security Act. 2001. **A, VI**
35. Advisory Committee on Dangerous Pathogens. Biological Agents: Managing the Risks in Laboratories and Healthcare Premises. Appendix 1.2 Transport of Infectious Substances - Revision. Health and Safety Executive 2008. **A, VI**
36. Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. MMWR Surveill Summ 2012;61:1-102. **B, V**
37. Health and Safety Executive. Control of Substances Hazardous to Health Regulations. The Control of Substances Hazardous to Health Regulations 2002 (as amended). Approved Code of Practice and guidance L5 (sixth edition). HSE Books,. 2013. **A, VI**
38. Health and Safety Executive. Five Steps to Risk Assessment: A Step by Step Guide to a Safer and Healthier Workplace. HSE Books,. 2002. **A, VI**
39. Health and Safety Executive. A Guide to Risk Assessment Requirements: Common Provisions in Health and Safety Law. HSE Books,. 2002. **A, VI**
40. Health Services Advisory Committee. Safe Working and the Prevention of Infection in Clinical Laboratories and Similar Facilities. HSE Books 2003. **A, VI**
41. British Standards Institution (BSI). BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets 2000. **A, VI**
42. British Standards Institution (BSI). BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets. Recommendations and guidance. 2005. 1-14. **A, VI**
43. Department of Health. Transport of Infectious Substances. Best Practice Guidance for Microbiology Laboratories. Department of Health. 1-13. 2007. **A, VI**
44. McLauchlin J. The pathogenicity of *Listeria monocytogenes*: a public health perspective. Rev Med Microbiol 1997;8:1-14.
45. Al-Tawfiq JA. *Listeria monocytogenes* bacteremia in a twin pregnancy with differential outcome: fetus papyraceus and a full-term delivery. JMicrobiolImmunolInfect 2008;41:433-6.
46. Rocourt J, Hof H, Schrettenbrunner A, Malinverni R, Bille J. [Acute purulent *Listeria seelingeri* meningitis in an immunocompetent adult]. SchweizMed Wochenschr 1986;116:248-51.
47. Lessing MP, Curtis GD, Bowler IC. *Listeria ivanovii* infection. JInfect 1994;29:230-1.
48. Cummins AJ, Fielding AK, McLauchlin J. *Listeria ivanovii* infection in a patient with AIDS. JInfect 1994;28:89-91.

49. Rocourt J, Schrettenbrunner A, Hof H, Espaze EP. [New species of the genus *Listeria*: *Listeria seeligeri*]. *PatholBiol(Paris)* 1987;35:1075-80.
50. Todeschini G, Friso S, Lombardi S, Casaril M, Fontana R, Corrocher R. A case of *Listeria murray/grayi* bacteremia in a patient with advanced Hodgkin's disease. *EurJClinMicrobiolInfectDis* 1998;17:808-10.
51. Rapose A, Lick SD, Ismail N. *Listeria grayi* bacteremia in a heart transplant recipient. *TransplInfectDis* 2008;10:434-6.
52. Brook I, Frazier EH. Significant recovery of nonsporulating anaerobic rods from clinical specimens. *Clin InfectDis* 1993;16:476-80.
53. McLauchlin J. The identification of *Listeria* species. *IntJFood Microbiol* 1997;38:77-81.
54. Barbuddhe SB, Maier T, Schwarz G, Kostrzewa M, Hof H, Domann E et al. Rapid identification and typing of *Listeria* species by matrix-assisted laser desorption ionization-time of flight mass spectrometry. *ApplEnvironMicrobiol* 2008;74:5402-7.
55. Hijazin M, Alber J, Lammler C, Weitzel T, Hassan AA, Timke M et al. Identification of *Trueperella* (*Arcanobacterium*) *bernardiae* by matrix-assisted laser desorption/ionization time-of-flight mass spectrometry analysis and by species-specific PCR. *JMed Microbiol* 2012;61:457-9.
56. Gasanov U, Hughes D, Hansbro PM. Methods for the isolation and identification of *Listeria* spp. and *Listeria monocytogenes*: a review. *FEMS MicrobiolRev* 2005;29:851-75.
57. Doumith M, Buchrieser C, Glaser P, Jacquet C, Martin P. Differentiation of the major *Listeria monocytogenes* serovars by multiplex PCR. *JClinMicrobiol* 2004;42:3819-22.
58. Halbedel S, Prager R, Fuchs S, Trost E, Werner G, Flieger A. Whole-Genome Sequencing of Recent *Listeria monocytogenes* Isolates from Germany Reveals Population Structure and Disease Clusters. *Journal of clinical microbiology* 2018;56.
59. Saltykova A, Wuyts V, Mattheus W, Bertrand S, Roosens NHC, Marchal K et al. Comparison of SNP-based subtyping workflows for bacterial isolates using WGS data, applied to *Salmonella enterica* serotype Typhimurium and serotype 1,4,[5],12:i. *PLOS ONE* 2018;13:e0192504.
60. Public Health England. Laboratory Reporting to Public Health England: A Guide for Diagnostic Laboratories. Public Health England 2016. 1-29. **A, VI**
61. Department of Health. Health Protection Legislation (England) Guidance. 1-112. 2010. **A, VI**
62. Scottish Government. Public Health (Scotland) Act. 2008. **A, VI**
63. Scottish Government. Public Health etc. (Scotland) Act 2008. Implementation of Part 2: Notifiable Diseases, Organisms and Health Risk States. 2009. **A, VI**
64. The Welsh Assembly Government. Health Protection Legislation (Wales) Guidance. 2010. **A, VI**
65. Home Office. Public Health Act (Northern Ireland) 1967 Chapter 36. 1967. **A, VI**