

COMMITTEE ON MUTAGENICITY OF CHEMICALS IN FOOD, CONSUMER PRODUCTS AND THE ENVIRONMENT (COM)

Follow up to Update of COM test guidance 3D Tissue Models for genotoxicity testing.

Background

1. The use of 3D models for genotoxicity testing has not previously been described in the COM (2011) “*Guidance On A Strategy For Genotoxicity Testing Of Chemical Substances*” (COM, 2011).
2. Development of 3D models is a rapidly evolving field, and therefore members considered it appropriate to prepare a document in this area, that can be updated at regular intervals, to support the Guidance Document. This discussion paper seeks to provide a brief summary of 3D models currently used for genotoxicity testing and those under development and/or validation.

Why 3D models are needed for genotoxicity testing

3. Although currently used *in vitro* genotoxicity testing batteries can reliably identify *in vivo* genotoxicants, there is a high degree of positives which, when tested *in vivo*, are non-genotoxic i.e. ‘false positives’. As a consequence, animal usage, testing time and costs are unnecessarily increased.
4. False positive findings are particularly associated with the *in vitro* mammalian monolayer cell assays such as the chromosomal aberration test, micronucleus test and mouse lymphoma assay (Kirkland et al., 2005; Kirkland et al., 2006; Kirkland et al., 2007; Matthews et al., 2006). This is considered to be due to a number of reasons, including the use of cell lines of rodent origin (V79, CHO or CL) that partially lack normal cell cycle control, have limited metabolic capacity (even with the addition of rat liver S9) and do not mimic site-specific metabolic capacity (Reus et al., 2013).
5. The main drivers for the development/use of 3D models were the Cosmetics Directive, which prevented the use of *in vivo* testing for cosmetics, and the 3Rs principle that requires the reduction in the use of animal toxicity testing. There are different types of 3D models, ranging from single cell microtissues to multi cell types grown within scaffolds. It is hoped the use of such models will reduce the number of false positives and improve the accuracy of predictions due to their improved metabolic capacity and the proximity to *in vivo* gene expression and protein function (Andres et al., 2012; Barcham et al., 2018).

3D models of skin

6. 3D models have, to date, mainly been developed for the skin. These models mimic the architectural features and behaviour of normal human skin and the changes that occur during early skin cancer progression and wound re-epithelialisation. Reconstructed 3D human epidermal skin models are used in OECD TG 431 (*in vitro* skin corrosion: Reconstituted human epidermis (RHE) test) (Kandárová et al., 2006; Kidd et al., 2007; OECD, 2016a) , which replaced the acute dermal irritation/corrosion test in rats (OECD TG 404). In addition, OECD TG 439 (*in vitro* skin irritation: Reconstituted human epidermis test) also utilises reconstructed 3D epidermal skin models (Alépée et al., 2010; Kandarova et al., 2009; OECD, 2015). Assessment of phototoxic properties (Jírová et al., 2005; Lelièvre et al., 2007) and sensitisation potential (dos Santos et al., 2011; Teunis et al., 2013) are also being explored using reconstructed 3D skin models and are considered to have a high potential to be accepted as OECD TGs (Reus et al., 2013).

7. For genotoxicity testing purposes, 3D skin models have been linked to the standard genotoxicity endpoints of the micronucleus test and Comet assay. Two endpoints are utilised to reflect different types of genetic damage, namely clastogenicity and aneugenicity and DNA strand breaks, incomplete repair sites and alkali labile sites, respectively. The 3D Skin Comet assay and Reconstituted Skin Micronucleus (RSMN) test are described in paragraphs 8 to 13 below. These assays allow the *in vitro* assessment of DNA damage following dermal exposure, which has only previously been possible using *in vivo* assays; this is despite dermal exposure being a common route for a number of compounds found in household products, cosmetics, and industrial chemicals (Reisinger et al., 2018).

3D Comet assay

8. The Comet assay has been adapted for use with two reconstructed full thickness human skin models: the EpiDerm™- and Phenion® FullThickness Skin Models. Both skin models are comprised of primary and p53 competent cells of human origin. These models have a number of advantages over current monolayer-type assays including: species specificity, with a phenotype close to native human skin; normal cell cycle control; DNA-repair competence; similar gene and protein expression patterns; and the mimicking of conditions of use for dermally applied substances/products (Reisinger et al., 2018).

9. As the Comet assay does not rely on proliferating cells and can be used with a wide range of cell types, it is particularly suitable for application to skin tissue models. The assay also detects a wide range of DNA damage including single-strand breaks from direct interaction of the test chemical or related to incomplete excision repair as well as alkali labile sites (OECD, 2016b). This ensures that both clastogenic DNA damage and lesions that may give rise to gene mutation are detected.

10. The 3D Skin Comet assay is currently undergoing inter-laboratory validation using the Phenion® Full-Thickness Skin Model to assess its potential use as a new *in vitro* tool for following up positive findings from the standard *in vitro* genotoxicity test battery for dermally applied chemicals. In the first phase using eight chemicals, the authors reported good intra- and inter-laboratory reproducibility with four of five laboratories reaching a 100% predictivity and the fifth yielding 70%. Further testing of 22 chemicals within a second validation phase is currently ongoing (Reisinger et al., 2018).

3D human reconstituted skin micronucleus assay

11. The RSMN assay has been developed to assess the genotoxicity of dermally applied compounds and utilises a highly differentiated *in vitro* model of the human epidermis (Episkin LM™) with automated micronucleus detection using the standard cytokinesis block micronucleus assay (Barcham et al., 2018). The RSMN offers a close approximation of natural human skin due to the origin of the cells used and its physiological properties. The model also allows topical administration which ensures that all parts of the model are exposed, regardless of the lipophilic nature of the test substance. In addition, the Episkin LM™ model has been shown to have a similar metabolic capacity to that of native human skin (Eilstein et al., 2014) allowing the assessment of genotoxic potential by metabolic activation as an intrinsic feature.

12. A global pre-validation of the assay is currently in progress (sponsored by the European Cosmetics Association and the European Center for Validation of Alternative Methods (ECVAM)). Initial findings suggest good inter-laboratory and inter-assay reproducibility for genotoxic chemicals that do not require metabolic activation (Aardema et al., 2013). However, further development of the assay is needed for genotoxic chemicals that require metabolism.

Other 3D tissue models

3D liver microtissue model

13. Conventional *in vitro* monolayer assays using hepatic cell lines may not be the most relevant assays to carry out functional and metabolic studies as the cells lose key liver specific functions, in particular cytochrome P450 activity (Godoy et al., 2013; Kim et al., 2011; Mingoia et al., 2007). In addition, non-parenchymal cells are absent which play an important role in clearance and in the initiation of an immune response. Due to the limited lifespan of the conventional assays, repeated exposures are not possible (Kermanizadeh et al., 2014).

14. A 3D liver microtissue model has been described (Messner et al., 2013; Kermanizadeh et al., 2014; Kratschmar DV, 2013) which has a number of advantages over conventional hepatic assays. These include: the use of primary human hepatic cells; viability of cells for long periods which allows multiple exposures to be assessed; maintenance of a high level of metabolic activity across the cells lifespan.

15. A 3D liver model utilising HepG2 cells grown using a 'hanging-drop' technique has been assessed for genotoxicity testing, with micronucleus detection in the 3D spheroid models. Micronucleus induction was seen to be greater in the 3D structures than in the 2D format (Shah et al., 2018).

3D tissue models of the airway epithelium

16. In conventional monolayer (2D) cultures of basal cells, only maintenance and expansion of cells is possible. However, in 3D airway tissue models, basal cells can differentiate into a mucociliary pseudostratified epithelium containing ciliated, goblet and basal cells. Other properties similar to the native human airway epithelium include beating cilia, mucus secretion, barrier properties and remodelling and restoration properties (Rock et al., 2009).

17. The two most widely used 3D tissue models of the airway epithelium are the air-liquid interface (ALI) cultures and sphere cultures.

18. ALI cultures reflect physiological conditions *in vivo*, with the respiratory epithelium being exposed to the air. These cultures are currently used to study cell biology and infection, culture patient-derived cells to model diseases, and test the effects of aerosolised particles (including drug formulations and cigarette smoke) on the respiratory epithelium (for example (Azzopardi et al., 2015)).

19. Sphere cultures are used in place of ALI cultures when a high throughput format is required. The process involves culture of airway epithelial cells in suspension which form spheres when placed into specialised coated tissue culture vessels. The cells differentiate to form a pseudostratified mucociliary epithelium with the apical surface of the cells pointing into the lumen (Rock et al., 2009).

Regulatory challenges

20. 3D human tissue models may offer an alternative testing strategy to *in vivo* assays for substances that are found to be positive using the traditional *in vitro* genotoxicity battery of tests.

21. Using historic data, chemicals that are positive for genotoxic activity *in vivo* have been shown to be positive in either the 3D-micronucleus or 3D-Comet assay. In the main, chemicals that are negative for genotoxic activity *in vivo* are also negative in the two 3D models (Kirkland et al., 2014).

22. However, even the most advanced of such models, the 3D Skin models, are currently only in the early stages of validation.

References

Aardema, M. J., Barnett, B. B., Mun, G. C., et al. 2013. Evaluation of chemicals requiring metabolic activation in the EpiDerm 3D human reconstructed skin micronucleus (RSMN) assay. *Mutat Res*, 750, 40-9.

- Alépée, N., Tornier, C., Robert, C., et al. 2010. A catch-up validation study on reconstructed human epidermis (SkinEthic™ RHE) for full replacement of the Draize skin irritation test. *Toxicology in Vitro*, 24, 257-266.
- Andres, E., Molinari, J., Remoue, N., et al. 2012. Successful micronucleus testing with the EPI/001 3D reconstructed epidermis model: preliminary findings. *Mutat Res*, 743, 36-41.
- Azzopardi, D., Haswell, L. E., Foss-Smith, G., et al. 2015. Evaluation of an air–liquid interface cell culture model for studies on the inflammatory and cytotoxic responses to tobacco smoke aerosols. *Toxicology in Vitro*, 29, 1720-1728.
- Barcham, R., Orsini, N., Andres, E., et al. 2018. Successful proof of concept of a micronucleus genotoxicity assay on reconstructed epidermis exhibiting intrinsic metabolic activity. *Mutat Res*, 829-830, 75-86.
- Com. 2011. *Guidance On A Strategy For Genotoxicity Testing Of Chemical Substances* [Online].
https://www.gov.uk/government/uploads/system/uploads/attachment_data/file/315793/testing_chemicals_for_genotoxicity.pdf. [Accessed].
- Dos Santos, G. G., Spiekstra, S. W., Sampat-Sardjoepersad, S. C., et al. 2011. A potential in vitro epidermal equivalent assay to determine sensitizer potency. *Toxicology in Vitro*, 25, 347-357.
- Eilstein, J., Lereaux, G., Budimir, N., et al. 2014. Comparison of xenobiotic metabolizing enzyme activities in ex vivo human skin and reconstructed human skin models from SkinEthic. *Arch Toxicol*, 88, 1681-1694.
- Godoy, P., Hewitt, N. J., Albrecht, U., et al. 2013. Recent advances in 2D and 3D in vitro systems using primary hepatocytes, alternative hepatocyte sources and non-parenchymal liver cells and their use in investigating mechanisms of hepatotoxicity, cell signaling and ADME. *Arch Toxicol*, 87, 1315-530.
- JíRová, D., Kejlová, K., Bendová, H., et al. 2005. Phototoxicity of bituminous tars—correspondence between results of 3T3 NRU PT, 3D skin model and experimental human data. *Toxicology in Vitro*, 19, 931-934.
- Kandarova, H., Hayden, P., Klausner, M., et al. 2009. In vitro skin irritation testing: Improving the sensitivity of the EpiDerm skin irritation test protocol. *Altern Lab Anim*, 37, 671-89.
- Kandárová, H., Liebsch, M., Spielmann, H., et al. 2006. Assessment of the human epidermis model SkinEthic RHE for in vitro skin corrosion testing of chemicals according to new OECD TG 431. *Toxicology in Vitro*, 20, 547-559.
- Kermanizadeh, A., Lohr, M., Roursgaard, M., et al. 2014. Hepatic toxicology following single and multiple exposure of engineered nanomaterials utilising a novel primary human 3D liver microtissue model. *Part Fibre Toxicol*, 11, 56.
- Kidd, D. A., Johnson, M. & Clements, J. 2007. Development of an in vitro corrosion/irritation prediction assay using the EpiDerm™ skin model. *Toxicology in Vitro*, 21, 1292-1297.
- Kim, B.-S., Park, I.-K., Hoshiba, T., et al. 2011. Design of artificial extracellular matrices for tissue engineering. *Progress in Polymer Science*, 36, 238-268.
- Kirkland, D., Aardema, M., Henderson, L., et al. 2005. Evaluation of the ability of a battery of three in vitro genotoxicity tests to discriminate rodent carcinogens and non-carcinogens I. Sensitivity, specificity and relative predictivity. *Mutat Res*, 584, 1-256.
- Kirkland, D., Aardema, M., Muller, L., et al. 2006. Evaluation of the ability of a battery of three in vitro genotoxicity tests to discriminate rodent carcinogens and non-

- carcinogens II. Further analysis of mammalian cell results, relative predictivity and tumour profiles. *Mutat Res*, 608, 29-42.
- Kirkland, D., Pfuhler, S., Tweats, D., et al. 2007. How to reduce false positive results when undertaking in vitro genotoxicity testing and thus avoid unnecessary follow-up animal tests: Report of an ECVAM Workshop. *Mutat Res*, 628, 31-55.
- Kirkland, D., Zeiger, E., Madia, F., et al. 2014. Can in vitro mammalian cell genotoxicity test results be used to complement positive results in the Ames test and help predict carcinogenic or in vivo genotoxic activity? I. Reports of individual databases presented at an EURL ECVAM Workshop. *Mutation Research/Genetic Toxicology and Environmental Mutagenesis*, 775-776, 55-68.
- Kratschmar Dv, M. S., Moritz W, Et Al. 2013. Characterization of a Rat Multi-Cell Type 3D-Liver Microtissue System. *J Tissue Sci Eng*, 4, 130 - 136.
- Lelièvre, D., Justine, P., Christiaens, F., et al. 2007. The episkin phototoxicity assay (EPA): Development of an in vitro tiered strategy using 17 reference chemicals to predict phototoxic potency. *Toxicology in Vitro*, 21, 977-995.
- Matthews, E. J., Kruhlak, N. L., Cimino, M. C., et al. 2006. An analysis of genetic toxicity, reproductive and developmental toxicity, and carcinogenicity data: I. Identification of carcinogens using surrogate endpoints. *Regulatory Toxicology and Pharmacology*, 44, 83-96.
- Messner, S., Agarkova, I., Moritz, W., et al. 2013. Multi-cell type human liver microtissues for hepatotoxicity testing. *Arch Toxicol*, 87, 209-13.
- Mingoia, R. T., Nabb, D. L., Yang, C. H., et al. 2007. Primary culture of rat hepatocytes in 96-well plates: effects of extracellular matrix configuration on cytochrome P450 enzyme activity and inducibility, and its application in in vitro cytotoxicity screening. *Toxicol In Vitro*, 21, 165-73.
- OECD 2015. *Test No. 439: In Vitro Skin Irritation: Reconstructed Human Epidermis Test Method*, Paris, France, OECD Publishing.
- OECD 2016a. *Test No. 431: In vitro skin corrosion: reconstructed human epidermis (RHE) test method. OECD Guidelines for the Testing of Chemicals, Section 4*. Paris: OECD Publishing.
- OECD 2016b. *Test No. 489: In Vivo Mammalian Alkaline Comet Assay*, Paris, France, OECD Publishing.
- Reisinger, K., Blatz, V., Brinkmann, J., et al. 2018. Validation of the 3D Skin Comet assay using full thickness skin models: Transferability and reproducibility. *Mutat Res*, 827, 27-41.
- Reus, A. A., Reisinger, K., Downs, T. R., et al. 2013. Comet assay in reconstructed 3D human epidermal skin models--investigation of intra- and inter-laboratory reproducibility with coded chemicals. *Mutagenesis*, 28, 709-20.
- Rock, J. R., Onaitis, M. W., Rawlins, E. L., et al. 2009. Basal cells as stem cells of the mouse trachea and human airway epithelium. *Proc Natl Acad Sci U S A*, 106, 12771-5.
- Shah, U. K., Mallia, J. O., Singh, N., et al. 2018. Reprint of: A three-dimensional in vitro HepG2 cells liver spheroid model for genotoxicity studies. *Mutat Res*, 834, 35-41.
- Teunis, M., Corsini, E., Smits, M., et al. 2013. Transfer of a two-tiered keratinocyte assay: IL-18 production by NCTC2544 to determine the skin sensitizing capacity and epidermal equivalent assay to determine sensitizer potency. *Toxicol In Vitro*, 27, 1135-50.

Glossary

3D Tissue Model	Artificially created environment in which biological cells are permitted to grow or interact with their surroundings in all three dimensions.
2D	Two dimensional
RSMN	Reconstituted Skin Micronucleus
ECVAM	European Center for Validation of Alternative Methods
HepG2 cells	Immortalised cell line consisting of human liver carcinoma cells

**NCET at WRc/IEH-C under contract supporting the PHE COT Secretariat
February 2019**