




Vernier Go Direct Charge Station	
Aspect	Details
Images	
PSD Number	2411-0064
Product Type	Electrical appliances and equipment - Custom Charging Dock
Product Identifiers	<p>Brand: Vernier Model/type: GDX-CRG (Also sold as part of a packaged item GDX-CRG, GDX-CRG-EX, GDX-CRG-NP, GDX-CRG-EURO, GDX-TMP-TP, GDX-TMP-TP-INTL, GDX-TMP-TP-NP, GDX-PH-TP, GDX-PH-TP-INTL, GDX-PH-TP-NP) Batch number: 17000001 - 23xxxxxx Order code: VR157352 (This product may also have been purchased as part of a Teacher Pack of sensors, order code VR157488 or VR157456. This recall applies to all Go Direct Charge Stations purchased prior to 15 December 2023.)</p>
Product Description	A charging station for Vernier wireless go direct sensors. Tools for use in educational laboratories – Schools / Colleges / Universities within the UK.
Country of Origin	Taiwan
Counterfeit	No
Risk Level	Not provided
Risk Type	Burns
Risk Description	If a product with an internal electrical short is plugged into this charge station the charge station could overheat to an unsafe temperature. The heat is high enough to melt the nearby ABS plastic of the GDX-CRG housing. Risk of burns. The product does not meet the requirements of the Electrical Equipment (Safety) Regulations 2016.
Corrective Measures	<p>The product has been recalled from end users.</p> <p>For further information, please access the Vernier recall alert here: https://www.vernier.com/support/gdx-crg-recall/?srsltid=AfmBOopq5_rfSDV-4IM0WPwfGMdwxs8oiHmKWbecQSBYTTCSSTITpfp.</p> <p>The importer has also been required to recall the products from end users. Please access the Instrument Direct Services recall site here https://www.indirectservices.co.uk/education/2024-vr157352-go-direct-charge-station-safety-recall/ and follow the steps to fill in the recall form.</p>
Online Marketplace	N/A
Notifier	Health and Safety Executive